



Acuo Admin Portal  
User Guide

# Copyright

Information in this document is subject to change without notice. The Acuo Admin Portal software (the "Software") described in this document is furnished only under a separate license agreement and may be used or copied only according to the terms of such agreement. It is against the law to copy the Software except as specifically allowed in the license agreement. This document or accompanying materials contains certain information which is confidential information of Hyland Software, Inc. and which is subject to the confidentiality provisions agreed to by you.

All data, names, and formats used in this document's examples are fictitious unless noted otherwise. Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright law, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the express written permission of Hyland Software, Inc.

©2018 Hyland Software, Inc. All rights reserved.

The Software may include software developed and copyrighted by third parties, including but not limited to the following:

Microsoft Corporation © Microsoft 2018

Nuance™ OCR © 2018 Nuance Communications, Inc. All rights reserved.

All rights reserved.

Further information regarding third-party software included in the Software can be found in the About box within the Software. Hyland, Hyland Software®, Acuo, and PACSgear are registered and/or unregistered trademarks of Hyland Software, Inc. in the United States and other countries. All other trademarks, service marks, trade names and products of other companies are the property of their respective owners.

# Table of Contents

<b>Copyright</b> .....	<b>2</b>
<b>What is Acuo Admin Portal?</b> .....	<b>7</b>
<b>Navigating Acuo Admin Portal</b> .....	<b>8</b>
View user information .....	8
Open a module .....	8
What is a system? .....	8
Select a system .....	8
Customize a list .....	8
<b>Patient Explorer</b> .....	<b>9</b>
What is Patient Explorer? .....	9
Find information by patient .....	9
Find information by study .....	9
Edit a patient record, study, or series .....	9
Delete information .....	10
Start quality control .....	10
Work with patient records .....	10
<i>Move a patient record</i> .....	10
<i>Merge patient records</i> .....	11
<i>Merge a study with a patient record</i> .....	11
Work with studies .....	11
<i>Move a study</i> .....	11
Work with series .....	12
<i>Move a series</i> .....	12
<i>Merge a series with a study</i> .....	12
Work with images .....	13
<i>Open a study in the NilRead viewer</i> .....	13
<i>Move an image</i> .....	13
<i>Merge an image with a series</i> .....	13
Search syntax .....	14
<i>Search syntax</i> .....	14
Table of search syntax .....	14
<b>Storage Management</b> .....	<b>15</b>
What is Storage Management? .....	15
What are archives, instances, and shares? .....	15

What is a storage application? .....	15
Add an archive device .....	15
Add a storage instance .....	16
Add a managed share .....	16
Add a storage application .....	16
Change the managed shares of a storage application .....	17
Edit or delete a configuration .....	17
Archive device settings .....	17
<i>Atmos</i> .....	17
<i>Azure</i> .....	19
<i>Bycast</i> .....	20
<i>Centera</i> .....	21
<i>DellCloud</i> .....	22
<i>DellDx</i> .....	23
<i>FileStore</i> .....	24
<i>Hitachi</i> .....	25
<i>Scality</i> .....	27
<b>Route Management</b> .....	<b>28</b>
What is Route Builder? .....	28
Work with routes .....	28
<i>What is a route?</i> .....	28
<i>Create a new called AE</i> .....	29
<i>Create an Acuo destination node</i> .....	29
<i>Create an external destination node</i> .....	29
<i>Create a route node</i> .....	30
<i>Create a complete route</i> .....	30
<i>Create a move route map</i> .....	30
<i>Edit route move and store options</i> .....	31
<i>View all routes of a destination</i> .....	31
Work with nodes .....	31
Edit node settings .....	31
Remove or delete a node .....	32
Work with destinations .....	32
<i>What is a destination?</i> .....	32
<i>Edit destination settings</i> .....	32
<i>Edit destination tuning parameters</i> .....	32

<b>Workflow</b>	<b>33</b>
What is Workflow?	33
View a list of workflow jobs	33
View a summary of store jobs	33
About performing actions on jobs	34
Perform an action on jobs	34
Available workflow job actions	34
Reschedule a move job	35
Change the destination of a store job	35
Resolve a store job conflict	35
<b>Enterprise Reporting</b>	<b>36</b>
What is Enterprise Reporting?	36
View and customize a report	36
Enterprise report customization options	36
<b>Worklist</b>	<b>40</b>
What is Worklist?	40
Edit or delete a patient	41
Patient settings	41
Edit or delete an event	42
Event settings	42
Add a worklist connection	44
Edit or delete a worklist connection	44
<b>XDS Explorer</b>	<b>45</b>
What is XDS Explorer?	45
Configure XDS	45
About PDQ and registry searches	45
Perform a PDQ search	45
Perform a registry search	46
Perform an MPQ search	46
Perform a VNA search	46
<b>System Management</b>	<b>46</b>
What is System Management?	46
Start a service	47
Stop a service	47
Restart a service	47
Configure Application Server with a load balancer	47

View a list of services .....	47
Restart all services .....	48
View system performance .....	48
<b>Lifecycle Management .....</b>	<b>48</b>
What is lifecycle management? .....	48
Create a retention policy .....	48
Schedule execution of a retention policy .....	49
<b>Auditing .....</b>	<b>49</b>
About Auditing .....	49
View an audit log .....	50
<b>System Configuration .....</b>	<b>50</b>
What is System Configuration? .....	50
Edit batch processing options .....	50
Batch processing options .....	50
Create an IPID domain .....	51
Edit or delete an IPID domain .....	51
Create an SOP .....	51
Edit or delete an SOP .....	52

# What is Acuo Admin Portal?

Acuo Admin Portal is a web-based program that provides administrators with tools to centrally manage medical information and systems across the enterprise. Acuo Admin Portal includes the following modules.

## **Patient Explorer**

View and edit patient demographics.

## **Storage Management**

Configure storage and retrieval of digital assets.

## **Route Management**

Create and configure routes between application entities (called AEs) and destinations.

## **Workflow**

Manage move and store jobs between systems.

## **Enterprise Reporting**

View analytics on VNA activity.

## **Worklist**

Manage HIS and RIS data.

## **XDS Explorer**

Browse XDS repositories.

## **System Management**

View, manage, or monitor application servers, Acuo Services and database performance.

## **Lifecycle Management**

Create data retention policies and schedule their execution.

## **Auditing**

View Acuo, IHE, user, and task audit logs.

## **System Configuration**

Manage and configure batch processing options, IPID domains and SOPs.

# Navigating Acuo Admin Portal

## View user information

Acuo Admin Portal offers a User Information page that allows you to view the following information: your name, username, the roles and group memberships to which you belong and the time you logged into Acuo Admin Portal. To view this information, complete one of the following steps.

- On the Acuo Admin Portal home screen, click your user name in top right corner of the screen, and select User Information. The system displays your User Information page.

## Open a module

Acuo Admin Portal comprises a set of modules and sub-modules in which you can perform a specific set of VNA functions. To open a module, complete one of the following steps.

1. On the Acuo Admin Portal home screen, click the button of the module that you want to open.
2. Optional. While in any module, in the upper left part of the screen, click the blue navigation arrow and select a module or sub-module.

## What is a system?

A system is a logical grouping of Acumen servers or an array of servers within Acuo Admin Portal that simplifies navigation and data management. A system contains one or more domains, and a domain contains one or more databases.

Systems are displayed in the Systems pane, but in certain contexts the Systems pane displays navigational items instead of systems.

## Select a system

In most modules, you can display a list of systems in the left pane. To select a system, complete the following steps.

1. To display the **Systems** pane, on the menu bar, select **Systems**.
2. Optional. In the **Search Systems** box, type any part of the system or domain name.
3. In the **Systems** pane, select a system.

## Customize a list

If a module displays a list such as a list of patients, jobs, or activity, you can often customize the list by sorting or applying filtering criteria. The customization options available differ depending on module and context. The following are commonly used customization options. For other options, see topics on lists in specific modules. To customize a list, complete the following steps.

1. On the menu bar, select from the available categories, enter filtering criteria, and then click **Apply**. For more information, see Search Syntax.
  - **Example.** In Patient Explorer click **Patient**, type a patient name, and then click **Apply**. You can enter criteria in more than one category.
2. Optional. On the column header, click a column name to sort by that column.



3. On the right side of the menu bar, select **Results per Page** and select a maximum number of results to display in the list.
4. On the right side of the menu bar, select **Columns** and then select a column to show or hide.

## Patient Explorer

### What is Patient Explorer?

Patient Explorer is a Acuo Admin Portal module that administrators can use to view and manage patient records and their associated studies, series, and images.

Common tasks in Patient Explorer include editing, removing or deleting, moving, and merging patient information. You can also view images at the study level. With Patient Explorer you can edit most DICOM metadata, but not raw image data.

### Find information by patient

You can enter one or more criteria to display matching patient records on the Patients tab. You can then view detailed information (such as demographics, studies, series, and images) in a patient record. To find information by patient, complete the following steps.

1. In **Patient Explorer**, select a system and select the **Patients** tab.
2. Optional. To customize the list of patients, see *Customize a list*.
3. To view patient studies, click the row containing the patient record.  
The studies appear under the patient row.
4. To view series within a study, click the row containing the study.
5. To view image metadata within a series, click the row containing the series.

### Find information by study

You can enter one or more criteria to display matching studies on the Studies tab. You can then view information (such as patient records, series, and images) in a particular study. To find information by study, complete the following steps.

1. In **Patient Explorer**, select a system and select the **Studies** tab.
2. Optional. To customize the list of studies, see *Customize a list*.
3. To view series within a study, click the row containing the study.  
The series appear under the study row.
4. To view image metadata within a series, click the row containing the series.

### Edit a patient record, study, or series

You can perform basic demographic editing tasks at the patient, study, or series level. For other quality control tasks, see "Quality Control." To edit a patient record, study, or series, complete the following steps.

1. In **Patient Explorer**, select a system.
2. On the **Patients** or **Studies** tab, find the patient, study, or series that you want to edit.

3. Click **Actions** and select **Edit**.
4. In the editing dialog box, make changes as necessary and click **Save**.

## Delete information

You can delete a patient record, study, series, or image from a system. For other quality control tasks, see "Quality Control." To delete information, complete the following steps.

1. In **Patient Explorer**, select a system.
2. On the **Patients** or **Studies** tab, find the patient record, study, series, or image that you want to delete.
3. Click **Actions** and select **Delete**.
4. In the confirmation dialog box, click **Delete [item]**.

## Start quality control

You can perform quality control tasks including the following.

- Modifying DICOM attributes
- Masking, shuttering, and labeling images
- Splitting series

For detailed instructions, see the Nil QC user's manual. To start quality control, complete the following steps.

1. In **Patient Explorer**, select a system.
2. On the **Patients** or **Studies** tab, find a patient, study, or series.
3. Click **Actions** and select **Start Quality Control**.

Nil QC launches and opens the associated patient record.

## Work with patient records

### Move a patient record

You can move (copy with a DICOM C-MOVE request) a patient record to one or more destinations through a previously configured route. To move a patient record, complete the following steps.

**Prerequisite** Configure at least one route in the Route Builder module, and a corresponding move route map. For more information, see [Create a move route map](#).

1. In **Patient Explorer**, select a system and select the **Patients** tab.
2. In the row of the patient that you want to move, click **Actions** and select **Move**.
3. In the **Move Patient** dialog box, in the **Move to Route** list, select a route.
4. Optional. Modify the **Priority** and **Scheduled Run Time** settings.
5. Click **Move Patient**.

## Merge patient records

You can merge one patient record (the source patient) with another patient record (the destination patient). This can be useful, for example, when a database contains more than one record for a patient, and you cannot delete either record because they contain non-identical information. To merge patient records, complete the following steps.

1. In **Patient Explorer**, select a system and select the **Patients** tab.
2. In the row containing the source patient, click **Actions** and select **Merge to Patient**.
3. To find a destination patient, complete the following substeps.
  1. In the **Patient Merge** dialog box, in the right pane, type a complete patient ID number in the **Patient ID** box, or type any part of the patient's first or last name in the **Patient Name** box, and then click **Search**.
  2. On the **Patients** tab, click the patient record that you want to be the destination patient.  
The name and ID of the patient record that you clicked appear in the left pane under Destination, and the **Studies** tab displays studies for the patient record. This is for reference only, and clicking a study, series, or image does not affect the patient merge operation.
4. Click **Merge Patient**.  
The patient records merge, and Acuo Admin Portal deletes the source patient.

## Merge a study with a patient record

You can merge a study from one patient record (the source study) with another patient record (the destination patient). To merge a study with a patient record, complete the following steps.

1. In **Patient Explorer**, select a system and select the **Studies** tab.
2. In the row containing the source study, click **Actions** and select **Merge to Patient**.
3. To find a destination patient record, complete the following substeps.
  1. In the **Study Merge** dialog box, type a complete patient ID in the **Patient ID** box, or type any part of the patient's first or last name in the **Patient Name** box, and then click **Search**.
  2. On the **Patients** tab, click the patient record that you want to be the destination patient.  
The name and ID of the patient record that you clicked appear in the left pane under Destination, and the **Studies** tab displays studies for the patient record. This is for reference only, and clicking a study, series, or image does not affect the study merge operation.
4. Click **Merge Study**.  
The source study merges with the destination patient, and Acuo Admin Portal deletes the source study.

## Work with studies

### Move a study

You can move (copy with a DICOM C-MOVE request) a study to one or more destinations through a previously configured route. To move a study, complete the following steps.

**Prerequisite** Configure at least one route by using the Route Builder module, and create a corresponding move route map. For more information, see [Create a move route map](#).

1. In **Patient Explorer**, select a system and select the **Studies** tab.
2. In the row containing the study to move, click **Actions** and select **Move**.
3. In the **Move Study** dialog box, in the **Move to Route** list, select a route.
4. Optional. Modify the **Priority** and **Scheduled Run Time** settings.
5. Click **Move Study**.

## Work with series

### Move a series

You can move (copy with a DICOM C-MOVE request) a series to one or more destinations through a previously configured route. To move a series, complete the following steps.

**Prerequisite** Configure at least one route by using the Route Builder module, and create a corresponding move route map. For more information, see [Create a move route map](#).

1. In **Patient Explorer**, select a system and select the **Studies** tab.
2. Find a series, click **Actions** and select **Move**.
3. In the **Move Series** dialog box, in the **Move to Route** list, select a route.
4. Optional. Modify the **Priority** and **Scheduled Run Time** settings.
5. Click **Move Series**.

### Merge a series with a study

You can merge a series from one study (the source series) with another study (the destination study). To merge a series with a study, complete the following steps.

1. In **Patient Explorer**, select a system and select the **Studies** tab.
2. Find a source series, click **Actions**, and then select **Merge to Study**.
3. To find a destination study, complete the following substeps.
  1. In the **Series Merge** dialog box, type a complete patient ID in the **Patient ID** box, or type any part of the patient's first or last name in the **Patient Name** box, and then click **Search**.
  2. On the **Patients** tab, click the patient record whose study you want to be the destination study.  
The name and ID of the patient record that you clicked appear in the left pane under Destination, and the **Studies** tab displays studies for the patient.
  3. On the **Studies** tab, click the study that you want to be the destination study.  
The UID of the study that you clicked appears in the left pane under Destination, and the **Series** tab displays series within the study. This is for reference only, and clicking a series or image does not affect the series merge operation.
4. Click **Merge Series**.  
The source series merges with the destination study, and Acuo Admin Portal deletes the source series.

## Work with images

### Open a study in the NilRead viewer

You can open a study in the NilRead viewer that comes with Acuo Admin Portal. To open a study in the viewer, complete the following steps.

1. In **Patient Explorer**, select a system and select the **Studies** tab.
2. In the row containing the study to open, click **Actions** and select **Launch in Viewer**.

### Move an image

You can move (copy with a DICOM C-MOVE request) an image to one or more destinations through a previously configured route. To move an image, complete the following steps.

**Prerequisite** Configure at least one route by using the Route Builder module, and create a corresponding move route map. For more information, see [Create a move route map](#).

1. In **Patient Explorer**, select a system and select the **Studies** tab.
2. Find a source image, click **Actions** and select **Move**.
3. In the **Move Instance** dialog box, in the **Move to Route** list, select a route.
4. Optional. Modify the **Priority** and **Scheduled Run Time** settings.
5. Click **Move Instance**.

### Merge an image with a series

You can merge an image from one series with another series. To merge an image with a series, complete the following steps.

1. In **Patient Explorer**, select a system and select the **Studies** tab.
2. Find a source image, click **Actions**, and then select **Merge to Series**.
3. To find a destination series, complete the following substeps.
  1. In the **Image Merge** dialog box, type a complete patient ID in the **Patient ID** box, or type any part of the patient's first or last name in the **Patient Name** box, and then click **Search**.
  2. On the **Patients** tab, click the patient record whose series you want to be the destination series. The name and ID of the patient record that you clicked appear in the left pane under Destination, and the **Studies** tab displays studies for the patient.
  3. On the **Studies** tab, click the study containing the series that you want to be the destination series. The UID (unique identifier) of the study that you clicked appears in the left pane under Destination, and the **Series** tab opens.
  4. On the **Series** tab, click the series that you want to be the destination series. The UID of the series that you clicked appears in the left pane under Destination, and the **Images** tab displays images within the series. This is for reference only, and clicking an image does not affect the image merge operation.

#### 4. Click **Merge Image**.

The source image merges with the destination study, and Acuo Admin Portal deletes the source image.

## Search syntax

### Search syntax

The following table describes the information that you can enter for each type of patient record or study search criterion.

Table of search syntax

Menu Bar	Item	Syntax
Patient	Patient ID	The exact patient ID.
	Issuer of Patient ID	Exact name of issuer.
	Patient Name	Any part of the patient's name. Case insensitive.
	Birth Date	YYYY-MM-DD. Enter beginning or ending date only, or both to specify a range.
Study	Accession Number	Exact accession number.
	Modality	Two-letter acronym, or multiple two-letter acronyms separated by backslashes. Case insensitive.
	Study UID	Exact study UID.
	Study Date	YYYY-MM-DD. Enter beginning or ending date only, or both to specify a range.

Searches also support the asterisk (\*) and question mark (?) wildcard characters.

- Gr?y matches Gray and Grey.
- 4\*8 matches 48, 408, 4018, etc.

# Storage Management

## What is Storage Management?

Storage Management is a Acuo Admin Portal module that administrators can use to configure storage and retrieval of digital assets. The module provides four areas of configuration: storage applications, managed shares, archive devices, and storage instances. You must configure these areas to enable subscriber clients such as AcuoMed to access the storage functionality of Acuo Admin Portal.

## What are archives, instances, and shares?

In the Storage Management module, you map your real-world storage systems to Acuo Admin Portal by configuring archive devices, storage instances, and managed shares. You must configure these before configuring a storage application.

### Archive device

An archive device is a mapping to a large-scale storage system.

### Managed share

A managed share is a mapping to a server or other machine-level storage space. You can assign multiple managed shares to a single machine.

### Storage instance

A storage instance is a mapping to a database instance on a managed share that stores images and other content. You can assign a storage instance to multiple managed shares.

## What is a storage application?

A storage application is a configuration that defines one storage instance, one archive device, or one of each. A storage application serves as the interface through which client applications such as AcuoMed access the storage functions of Acuo Admin Portal.

When you add a storage application for a managed share, you assign the managed share by selecting a storage instance. You must add a storage application for each storage instance or archive that you want to be available to clients.

Storage applications are also available for use in other modules of Acuo Admin Portal. For example, in Route Builder you can create an Acuo Destination node that points to a storage application.

## Add an archive device

You must add an archive device for each real-world archive that you want to be available in Acuo Admin Portal. To add an archive device, complete the following steps.

1. In **Storage Management**, in the left pane, select **Archive Devices**.
2. On the **Archive Devices** page, click **Add New Archive Device** and select a system in the list.
3. In the dialog box that appears, configure the archive device settings.

**Note:** For a description of the available settings, see Archive device settings.

4. Click **Add New Archive Device**.

## Add a storage instance

A storage instance identifies a database and its host machine to Acuo Admin Portal for purposes of image storage. You define the path to the storage instance when you add a managed share. To add a storage instance, complete the following steps.

1. In **Storage Management**, in the left pane, click **Storage Instances**.
2. On the **Storage Instances** page, click **Add Extended Storage Instance**.
3. In the **Add Extended Storage Instance** dialog box, in the **Data Source** box, type the name of the machine to host the database instance.
4. Under **Database Name**, do one of the following:
  - Type the name of a database that exists on the data source.
  - To create a new database, type a name for the database and select the **Create a New Database** check box.
5. Click **Add Extended Storage Instance**.

## Add a managed share

You must add a managed share for each storage area on a machine that you want to use in Acuo Admin Portal. To add a managed share, complete the following steps.

1. In **Storage Management**, in the left pane, select **Managed Shares**.
2. On the **Managed Shares** page, click **Add New Managed Share**.
3. In the **Add New Managed Share** dialog box, complete the following substeps.
  1. In the **Name** box, type a name to identify the managed share in Acuo Admin Portal.
  2. In the **Path** box, type the path to the storage area in the format `\\[server name][directory]`.
  3. In the **Storage Instance** list, select a storage instance for the managed share.

**Note:** Storage Management locates the storage instance in the path that you specify.

4. Click **Add New Managed Share**.

## Add a storage application

**Prerequisite** Add a storage instance and managed share, or an archive.

You must add a storage application to make managed shares and archives available to Acuo Admin Portal clients. To add a storage application, complete the following steps.

1. In **Storage Management**, in the left pane, select **Storage Applications**.
2. On the **Storage Applications** page, click **Add New Storage Application**.
3. In the **Add New Storage Application** dialog box, complete the following substeps.
  1. In the **Application Name** box, type a name to identify the storage application in Acuo Admin



Portal.

2. Under **Type**, click the button that corresponds to the type of data that the application stores.

**Note:** Generic refers to non-DICOM images.

3. In the **Storage Instance** and **Archive** lists, select a storage instance, archive, or both.
4. Click **Add New Storage Application**.

## Change the managed shares of a storage application

You can assign additional managed shares to an existing storage application, or remove its assigned shares. To change the managed shares of a storage application, complete the following steps.

1. In **Storage Management**, in the left pane, select **Storage Applications**.
2. In the right pane, in the row containing a storage application, click **Shares**.
3. In the **Shares** dialog box, in the list of available shares, select a managed share to switch its status between unassigned and assigned.

**Note:** Assigned shares appear in blue with a white check mark.

4. In the **Active Share** list, select the managed share in which you want to store new files.
5. Click **Update Shares**.

## Edit or delete a configuration

You can edit or delete the configuration of an archive device, storage instance, managed share, or storage application. To edit or delete a configuration, complete the following steps.

1. In **Storage Management**, in the left pane, select a system.
2. In the right pane, in the row containing the configuration that you want to edit or delete, do one of the following:
  - To delete the configuration, click **Delete**.
  - To edit the configuration, click **Edit**, edit the configuration as needed, and then click the update button.

## Archive device settings

The following tables describe the settings, by vendor, that are available when you add an archive device. See "Add an archive device" for the configuration procedure.

### Atmos

The following settings are available for an Atmos archive device.

<b>Property Type</b>	<b>Property Value</b>	<b>Property Description</b>
Archive Device Settings	Name	A name to identify the archive device in Storage Management.
	Connection String	The IP address of the archive device
	Port	The port number on which Acuo Admin Portal can communicate with the archive device.
	Use SSL	Identifies whether archive device uses SSL encryption.
Security Settings	Username and Password	The login credentials for the archive device.
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer Launch Time	The time at which Storage Management starts storing data to the archive device. Set this time so that it does not overlap with your SQL maintenance job schedule.
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.

Property Type	Property Value	Property Description
	Maximum Addresses	The maximum number of addresses.

## Azure

The following settings are available for an Azure archive device.

Property Type	Property Value	Property Description
Archive Device Settings	Name	A name to identify the archive device in Storage Management.
	Connection String	The IP address of the archive device
	Port	The port number on which Acuo Admin Portal can communicate with archive device.
	Use SSL	Identifies whether the archive device uses SSL encryption.
Security Settings	Account Name and Account Key	The login credentials for the archive device.
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer Launch Time	The time at which Storage Management starts storing data to the archive device. Set this time so that it does not overlap with your SQL maintenance job schedule.
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.

Property Type	Property Value	Property Description
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.
	Maximum Addresses	The maximum number of addresses.

## Bycast

The following settings are available for a Bycast archive device.

Property Type	Property Value	Property Description
Archive Device Settings	Name	A name to identify the archive device in Storage Management.
	Connection String	The IP address of the archive device
	Port	The port number on which Acuo Admin Portal can communicate with the archive device.
	Query Retrieve Port	The port number on which Acuo Admin Portal can receive queries and retrieve results from the archive device.
	Use SSL	Identifies whether the archive device uses SSL encryption.
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer Launch Time	The time at which Storage Management starts storing data to the archive device. Set this time so that it does not overlap with your SQL maintenance job schedule.

Property Type	Property Value	Property Description
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.
	Maximum Addresses	The maximum number of addresses.

## Centera

The following settings are available for a Centera archive device.

Property Type	Property Value	Property Description
Archive Device Settings	Name	A name to identify the archive device in Storage Management.
	Connection String	The IP address of the archive device
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer	The time at which Storage Management starts storing data to the archive

Property Type	Property Value	Property Description
	Launch Time	device. Set this time so that it does not overlap with your SQL maintenance job schedule.
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.
	Maximum Addresses	The maximum number of addresses.

## DellCloud

The following settings are available for a DellCloud archive device.

Property Type	Property Value	Property Description
Archive Device Settings	Name	A name to identify the archive device in Storage Management.
	Connection String	The IP address of the archive device
	Port	The port number on which Acuo Admin Portal can communicate with archive device.

Property Type	Property Value	Property Description
	Use SSL	Identifies whether the archive device uses SSL encryption.
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer Launch Time	The time at which Storage Management starts storing data to the archive device. Set this time so that it does not overlap with your SQL maintenance job schedule.
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.
	Maximum Addresses	The maximum number of addresses.

## DelIDx

The following settings are available for a DelIDx archive device.

Property Type	Property Value	Property Description
Archive Device	Name	A name to identify the archive device in Storage Management.

Property Type	Property Value	Property Description
Settings		
	Connection String	The IP address of the archive device
	Port	The port number on which Acuo Admin Portal can communicate with archive device.
	Use SSL	Identifies whether the archive device uses SSL encryption.
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer Launch Time	The time at which Storage Management starts storing data to the archive device. Set this time so that it does not overlap with your SQL maintenance job schedule.
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.
	Maximum Addresses	The maximum number of addresses.

## FileStore

The following settings are available for a FileStore archive device.



<b>Property Type</b>	<b>Property Value</b>	<b>Property Description</b>
Archive Device Settings	Name	A name to identify the archive device in Storage Management.
	NAS Node Name	The name of the NAS node that Storage Management can use.
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer Launch Time	The time at which Storage Management starts storing data to the archive device. Set this time so that it does not overlap with your SQL maintenance job schedule.
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.
	Maximum Addresses	The maximum number of addresses.

## Hitachi

The following settings are available for a Hitachi archive device.

<b>Property Type</b>	<b>Property Value</b>	<b>Property Description</b>
Archive Device Settings	Name	A name to identify the archive device in Storage Management.
	Connection String	The IP address of the archive device
	Port	The port number on which Acuo Admin Portal can communicate with archive device.
	Use SSL	Identifies whether the archive device uses SSL encryption.
Security Settings	Username and Password	The login credentials for the archive device.
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer Launch Time	The time at which Storage Management starts storing data to the archive device. Set this time so that it does not overlap with your SQL maintenance job schedule.
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.

Property Type	Property Value	Property Description
	Maximum Addresses	The maximum number of addresses.

## Scality

The following settings are available for a Scality archive device.

Property Type	Property Value	Property Description
Archive Device Settings	Name	A name to identify the archive device in Storage Management.
	Connection String	The IP address of the archive device
	Port	The port number on which Acuo Admin Portal can communicate with archive device.
	Use SSL	Identifies whether the archive device uses SSL encryption.
Security Settings	Account Name and Account Key	The login credentials for the archive device.
	Bucket Name	The name of the bucket that is indexed for file storage.
Archive Writer Options	Writer Interval	The time interval during which Storage Management continuously writes files to the archive device.
	Writer Launch Time	The time at which Storage Management starts storing data to the archive device. Set this time so that it does not overlap with your SQL maintenance job schedule.
	Maximum Threads	The total number of threads to the archive device that can be open simultaneously.
Tuning Parameters	Days Old	The number of days to wait before storing files to the archive device.

Property Type	Property Value	Property Description
	Days on Cache	The number of days after which Storage Management removes files from the cache. This setting must be greater than the Days Old setting.
	Maximum Clip Size	The maximum container size, in megabytes.
Advanced Options	Writer Retry Count	The number of times Storage Management attempts to write data.
	Writer Retry (Minutes)	The length of time between writer retries.
	Maximum Addresses	The maximum number of addresses.

## Route Management

### What is Route Builder?

Route Builder is a sub-module of the Route Manager Acuo Admin Portal module that administrators can use to configure routes for sending data from a called AE to one or more destinations.

Route Builder graphically depicts routes as a series of connected nodes representing called AEs, route names, and destinations. Displaying routes in this way makes it easier to understand the flow of information in a route, and to modify settings for each of its devices.

### Work with routes

#### What is a route?

When you perform a move operation in Patient Explorer, information flows along a route from a called AE to one or more destinations. Route Builder displays the route as a series of connected nodes representing the called AE, the route name, and the destinations.

There are four types of routes.

#### Basic route

Also referred to as a primary route, this is a standard route between an AE and one or more destinations. Move jobs on a basic route remain in a queue until reconciliation, so there may be a delay until the jobs arrive at the destinations.

#### Stat route

An expedited route that bypasses reconciliation so that jobs arrive at their destinations immediately.

#### Tag failure route

A conditional route that move jobs take when they satisfy user-defined rules.

### **Auto Change Dest.**

A special type of route in which Route Builder can re-route store jobs in the case of a UID failure. You can only create one route of this type.

## Create a new called AE

You can create a new called AE for use in creating routes. To create a new called AE, complete the following steps.

1. In **Route Builder**, on the menu bar, select **Create New Called AE**.
2. In the **Create New Called AE** dialog box, complete the following substeps.
  1. In the **Called AE Title** box, type the AE title of the device.
  2. In the **Local SCP** list, select an IP address and port number of a service class provider.
  3. Click **Create New Called AE**.

## Create an Acuo destination node

You can create an Acuo destination node for use in creating routes. For more information, see [What is a destination?](#) To create an Acuo destination node, complete the following steps.

1. In **Route Builder**, on the menu bar, select **Showing**, select the **Routes** tab, and then select a route on which to create the destination.
2. In the route display area, click the route to which you want to add the destination.
3. In the settings area, under **Acuo Destinations**, click **Create a New Destination**.
4. In the **Create New Acuo Destination** dialog box, type a name for the destination, and then do one of the following:
  - To use an existing database for the destination, click **Choose an Existing DicomDb**.
  - To configure a new database for the destination, click **Create a New DicomDb**, and then enter settings for the new database.
5. In the **DicomDb** list, select a database.
6. In the **Acuo Store Application Name** list, select your Acuo Store application, and then click **Add Destination**.

## Create an external destination node

You can create an external destination node for use in creating routes. For more information, see [What is a destination?](#) To create an external destination node, complete the following steps.

1. In **Route Builder**, on the menu bar, select **Showing**, select the **Routes** tab, and then select a route on which to create the destination node.
2. In the route display area, click the route to which you want to add the destination.
3. In the settings area, under **External Destinations**, click **Create a New Destination**.
4. In the **Create New External Destination** dialog box, complete the following substeps.

1. In the **Destination Name** box, type a name.
2. In the **Called AE** and **Calling AE** boxes, type the called AE title and calling AE title of the destination device.
3. In the **Host** and **Port** boxes, type the host name and port number of the destination device.
4. Click **Create New External Destination**.

## Create a route node

You can create a route node for use in creating complete routes. To create a route node, complete the following steps.

1. In **Route Builder**, on the menu bar, select **Create New Route**.
2. In the **Create New Route** dialog box, type a name for the route node and select the route type.
3. Click **Create New Route**.

**Note:** The route node appears in the route display area. You can add destination nodes to the route node here, but to create a complete route, see [Create a complete route](#).

## Create a complete route

The main purpose of Route Builder is to create a complete route that includes nodes for called AEs, routes, and destinations. To create a complete route, complete the following steps.

**Prerequisite** Create at least one new called AE node, one route node, and one destination node.

1. In **Route Builder**, on the menu bar, select **Showing**, select the **Called AEs** tab, and then select an AE to be the first node in the complete route.
2. Under **Primary Route**, **Stat Route**, or **Tag Routing**, click **Assign an Existing Route**.
3. In the **Select Route** dialog box, select a route, and then click **Select Route**.
4. In the route display area, select the route node.
5. In the settings area, click **Add Existing Destinations**.
6. In the **Add Destinations** dialog box, select one or more destinations and then click **Add Destinations**.

## Create a move route map

By default, routes you create in Route Builder are store routes. If you want to move records on a route, you must create a move route map. To create a move route map, complete the following steps.

1. In **Route Builder**, in the **Systems** pane, select **Move Route Maps**.
2. On the menu bar, select **Create Move Route Map**.
3. In the **Create Move Route Map** dialog box, complete the following substeps.
  1. In the **Map Name** box, type a name for the move route map.
  2. In the **Destination AE** box, type the destination for move jobs.

**Note:** This corresponds to the DICOM Move AE.

3. In the **Route** list, select the route that you want to be the move route.

#### 4. Click **Create Move Route Map**.

### Edit route move and store options

Administrators can edit advanced move and store options for a route. To edit route move and store options, complete the following steps.

1. In **Route Builder**, in the main viewing area, click a route node.
2. In the settings area, click **Move and Store Options**.
3. In the **Route Options** dialog box, complete the following substeps.
  1. Under **Move Level**, select whether AcuoMed performs move operations at the patient, study, series, or image level.
  2. To prevent unnecessary move operations, select the **Only move assets not found at target destinations** check box.
  3. To issue a C-FIND prior to a move operation if the move levels of the source and route are the same, select the **Query sources when move and route levels are the same** check box.
  4. To prevent issuing of C-STORE operations, such as when an archive is present, and to bring offline images online, select the **Do not issue C-STOREs from Local DBs** check box.
  5. To limit the number of responses to a find request, in the **Find Response Limit** box, type a number.
  6. Click **Save**.

### View all routes of a destination

To anticipate the system-wide impact that changing a destination might have, you can display all routes that include a specific destination. To view all routes of a destination, complete the following step.

1. In **Route Builder**, on the menu bar, select **Showing**.
2. Select the **Destinations** tab, and then select a destination.

**Result** All routes containing the destination appear in the route display area.

### Work with nodes

#### Edit node settings

You can edit settings for any called AE, route, or destination node on a route. To edit node settings, complete the following steps.

1. In **Route Builder**, on the route display area, click a node.
2. In the settings area, click **Edit [node] Settings**.
3. In the settings dialog box for the node, edit settings and then click **Save**.

## Remove or delete a node

Removing a node removes it from the current route. Deleting a node deletes it from all routes, and makes it unavailable for future use. You can only delete, and not remove, a called AE node. To remove or delete a node, complete the following steps.

1. In **Route Builder**, in the route display area, click a node.
2. In the settings area, do one of the following:
  - To remove the node, select **Remove [Node Type] from [Parent Node]**.
  - To delete a node, select **Delete [Node Type]**.

## Work with destinations

### What is a destination?

A destination is a DICOM device that you assign to a route for receiving information from a called AE.

You can configure two types of destinations.

#### Acuo destination

An AE on the local network, or a local DICOM destination that you previously configured on the current system.

#### External destination

A non-local DICOM device such as another Acumen system, PACS, VNA, viewer, archive, or laptop.

### Edit destination settings

You can edit the settings that you configured when creating a destination node. To edit destination settings, complete the following steps.

1. In **Route Builder**, in the route display area, click the destination node whose settings you want to edit.
2. In the settings area, select **Edit Destination Settings**.
3. In the dialog box that appears, edit settings and click **Save**.

### Edit destination tuning parameters

You can edit destination tuning parameters to adjust the number of database commands available for DICOM C-STORE and C-FIND operations. To edit tuning parameters, complete the following steps.

1. In **Route Builder**, in the route display area, click the destination whose tuning parameters you want to edit.
2. In the settings area, select **Tuning Parameters**.
3. In the tuning parameters dialog box, under **C-STORE Options** or **C-FIND Options**, complete the following substeps.
  1. In the **Maximum Connections** box, type the maximum number of commands that the database can receive.



2. In the **Minimum Free Connections** box, type a number. When the number of free connections falls to this number, the system creates new free connections, up to the maximum.
3. In the **Maximum Free Connections** box, type the maximum number of free connections that the system can create.
4. Click **Save**.

## Workflow

### What is Workflow?

Workflow Management (Workflow) is a Acuo Admin Portal module that administrators can use to manage move and store jobs that are queued by other Acuo Admin Portal modules. Workflow contains two sub-modules, Move Manager and Store Manager.

#### Move Manager

Displays lists of active and completed move jobs, and enables you to pause, restart, or reschedule them.

#### Store Manager

Displays aggregate statistics on active and completed store jobs. It also displays system-specific lists that provide various job management functions including prioritization, hold, pause, cancel, destination change, and conflict resolution.

### View a list of workflow jobs

You can monitor workflow activity by displaying lists of workflow jobs.

1. To view a list of workflow jobs, do one of the following:
  - To view move jobs, in **Move Manager**, select a system.
  - To view store jobs, in **Store Manager**, select a system other than a top-level system (such as the default system).
2. Optional. Customize the list. For more information, see [Customize a list](#).

### View a summary of store jobs

You can view a summary of active and completed store jobs that includes status, average queue and run times, and exceptions and conflicts. To view a summary of store jobs, complete the following steps.

1. In **Store Manager**, select a top-level system such as the default system.

**Note:** On the summary screen, you can move the pointer over graphs to view details.

2. Optional. Do one or more of the following:
  - To view a list of jobs by status or other attributes, in the **Active Jobs** or **Recently Completed Patients** header, click a number under a column.
  - To view a list of jobs by error type, under **Action Center**, select an error type.

## About performing actions on jobs

When you select a system in a Workflow module, the name of the actions menu on the menu bar changes accordingly. For example, depending on which system you select the menu could be Queue Actions, Destination Actions, or Patient Actions. Additionally, the actions available in the menu change depending on the selected system and the types of job currently displayed.

Depending on context, workflow actions available for multiple jobs may also be available for individual jobs on their Actions menus.

## Perform an action on jobs

You can perform actions on jobs such as changing their status, schedule, or destinations. To perform an action on jobs, complete the following steps.

1. In **Move Manager** or **Store Manager**, select a system.
2. Do one of the following:
  - To modify multiple jobs, on the menu bar, select the action menu displayed and select an action.
  - To modify a single job, in the job row, click **Action** and select an action.

**Note:** For a description of available actions, see Available workflow job actions.

## Available workflow job actions

The following list describes the actions that are available in Move Manager, Store Manager, or both.

### **Pause**

Pauses one or more jobs. Acuo Admin Portal also pauses jobs automatically when they experience a conflict.

### **Expedite**

Changes the status of jobs to Expedite. Acuo Admin Portal processes these jobs first.

### **Lower Priority**

Changes the status of jobs to Lower Priority. Acuo Admin Portal processes these jobs last.

### **Resume**

Resumes jobs that are paused or on hold.

### **Hold**

Places paused jobs on hold for further processing.

### **Cancel**

Cancels jobs that are paused or on hold and removes them from the queue.

### **Reschedule**

Reschedule move jobs to change their time of execution. For more information, see [Reschedule a move job](#).

### Change Destination

Changes the destinations that you pre-configured in Route Builder. See [Change the destination of a store job](#).

### Resolve Conflict

Resolves conflicts in paused jobs. See [Resolve a store job conflict](#).

## Reschedule a move job

You can change the start date and time of move jobs that have not yet started. To reschedule a move job, complete the following steps.

1. In **Move Manager**, select a system.
2. In the list of jobs, find the job to reschedule, click **Actions**, and then select **Reschedule**.
3. In the **Reschedule Job** dialog box, select a new date and time and click **Reschedule Job**.

## Change the destination of a store job

A store job moves along a predetermined route, but you can change its destinations while the job is still active. To change the destination of a store job, complete the following steps.

1. In **Store Manager**, select a system.
2. In the row containing the job to modify, click **Actions** and select **Change Destination**.
3. In the **Change Destination** dialog box, in the **New Destination** list, select a new destination and click **Change Destination**.

## Resolve a store job conflict

**Prerequisite** AcuoMed 6.0.3 HF1

During store operations, conflicts can arise when patient or study information (such as patient IDs or study UUIDs) differs between the source and store destination. Store manager detects conflicts and provides the user with options for resolving them. To resolve a store job conflict, complete the following steps.

1. In **Store Manager**, select a system.
2. In the list of store jobs, find a job with a conflict.

**Note:** Store Manager displays a red icon next to jobs with conflicts.

3. Click **Actions** and select **Resolve Conflict**.
4. On the conflict summary page, click **Next**.
5. On the **Choose the Appropriate Conflict Resolution** page, select a conflict resolution option and click **Next**.

6. On the **Confirm Conflict Resolution** page, click **Next**.

Store Manager executes the conflict resolution option and displays the result.

## Enterprise Reporting

### What is Enterprise Reporting?

Enterprise Reporting is a Acuo Admin Portal module that administrators can use to view lists, charts, and graphs of VNA activity. It comprises the following sub-modules.

#### **Route Graph**

Displays a graph of routes on multiple systems to help you visualize your overall routing scheme.

#### **Volume Analysis**

Displays color-coded charts and graphs to help you visualize the volume of traffic by route and modality.

#### **Standard Reports**

Displays lists of workflow and archiving activity customizable by destination, AE, server, and database.

#### **Archiving Analysis**

Displays system archive reports in chart and list format.

#### **Workflow Analysis**

Displays workflow activity by AE, modality, destination, and server in chart and list format.

#### **Storage Conflicts**

Displays a graph of storage conflicts by AE, conflict type, destination, and server.

#### **Study Usage**

Displays a graph of study usage over time.

### View and customize a report

To view and customize an enterprise report, complete the following steps.

1. In **Enterprise Reporting**, click a sub-module and select a system.
2. To customize the report, use the available menu bar options (see "Enterprise report customization options" and "Customize a list").

### Enterprise report customization options

The following tables describe the customization options that are available in each Enterprise Reporting module.

#### **Route Graph**

Tool	Option	Description
Filter	All Routes	Displays all routes
	Complete Routes	Displays routes with at least 3 nodes
	Incomplete Routes	Displays routes with fewer than 3 nodes

### Volume Analysis

Tool	Option	Description
Date Range	<ul style="list-style-type: none"> <li>• Begin/End</li> <li>• Last 7/30/90 days</li> </ul>	Displays only activity occurring within the specified date range
Previous Month/Week		Displays data one month or week before the data currently displayed
Next Month/Week		Displays data one month or week after the data currently displayed

### Standard Reports

Tool	Option	Description
Date Range	<ul style="list-style-type: none"> <li>• Begin/End</li> <li>• Last 7/30/90 days</li> </ul>	Displays only activity occurring within the specified date range
Report	Jobs Summary	Displays a list of jobs
	Destination Volume	Displays number of jobs by destination
	AE Volume	Displays number of jobs by AE
	Server Volume	Displays number of jobs by server
	Database Volume	Displays job detail by server and database
Daily/Monthly/Yearly		Displays jobs by server since the previous day, month, or year. Clicking these automatically modifies the date range setting.

<b>Tool</b>	<b>Option</b>	<b>Description</b>
Date	Context dependent	Filters the current report by date
Destinations	Context dependent	Filters the current report by destination
Server Name	Context dependent	Filters the current report by server
Download to CSV		Downloads the current report as a CSV file
Print		Prints the current report

### Archiving Analysis

<b>Tool</b>	<b>Option</b>	<b>Description</b>
Date Range	<ul style="list-style-type: none"> <li>• Begin/End</li> <li>• Last 7/30/90 days</li> </ul>	Displays only activity occurring within the specified date range
Patient		Sorts the current report by patient
Study		Sorts the current report by study
Series		Sorts the current report by series
Image		Sorts the current report by image
Size		Sorts the current report by size
Daily/Monthly/Yearly		Displays jobs by server since the previous day, month, or year. Clicking these automatically modifies the date range setting.
Chart		Displays the report in chart form
Size/Byte		Specifies the size you want the output to display.
Data		Displays the report as a list
Date	Context dependent	Filters the current report by date
DicomDb	Context dependent	Filters the current report by DICOM database

<b>Tool</b>	<b>Option</b>	<b>Description</b>
AcuoStore	Context dependent	Filters the current report by AcuoStore instance
Download to CSV		Downloads the current report as a CSV file
Print		Prints the current report

### **Workflow Analysis**

<b>Tool</b>	<b>Option</b>	<b>Description</b>
Date Range	<ul style="list-style-type: none"> <li>• Begin/End</li> <li>• Last 7/30/90 days</li> </ul>	Displays only activity occurring within the specified date range
AEs		Displays the number of move/store jobs by AE
Modalities		Displays the number of move/store jobs by modality
Destinations		Displays the number of move/store jobs by destination
Servers		Displays the number of move/store jobs by server
Daily/Monthly/Yearly		Displays jobs by server since the previous day, month, or year. Clicking these automatically modifies the date range setting.
Chart		Displays the report in chart form
Data		Displays the report as a list
Date	Context dependent	Filters the current report by date
Modality	Context dependent	Filters the current report by modality
Download to CSV		Downloads the current report as a CSV file
Print		Prints the current report

## Workflow Analysis

Tool	Option	Description
Date Range	<ul style="list-style-type: none"><li>• Begin/End</li><li>• Last 7/30/90 days</li></ul>	Displays only activity occurring within the specified date range
AEs		Displays the number of conflicts by AE
Conflicts		Displays the number of conflicts by conflict type
Destinations		Displays the number of conflicts by destination
Servers		Displays the number of conflicts by server
Daily/Monthly		Displays jobs by server since the previous day, month, or year. Clicking these automatically modifies the date range setting.

## Study Usage

Tool	Option	Description
Date Range	<ul style="list-style-type: none"><li>• Begin/End</li><li>• Last 7/30/90 days</li></ul>	Displays only activity occurring within the specified date range
Download to CSV		Downloads the current report as a CSV file

## Worklist

### What is Worklist?

Worklist is a Acuo Admin Portal module that administrators can use to manage scheduled patients, worklist events and worklist configuration and prefetching rules. Worklist contains three sub-modules, Patients, Events and Configuration.

#### Patients

Displays a list of scheduled patients, and enables you to view, delete, and edit patients.

#### Events

Displays a list of patient events, and enables you to view, delete, and edit events.

#### Configuration



Displays a list of worklist connections, and enables you to create and modify new worklist connections.

## Edit or delete a patient

To edit or delete a patient, complete the following steps.

1. In **Worklist**, in the left pane, click **Patients**.
2. In the right pane, in the row containing the patient that you want to edit or delete, do one of the following:
  - To delete a patient, click **Delete** and in the **Delete Patient** dialog box, click **Delete Patient** to confirm. Important: Deleting a patient also deletes all associated events.
  - To edit a patient, click the patient's name and in the **Edit Patient** dialog box, edit the event as needed, and then click **Save Changes**. For a description of the available settings, see "Patient Settings".

## Patient settings

The following table describes the settings that are available when you edit a patient. See "Edit or delete a patient" for the configuration procedure.

Property Type	Property Value	Property Description
Patient Information	Patient Name	The patient's full name.
	Patient ID	The unique identifier for the patient.
	Issuer of Patient ID	The unique identifier of the individual that issued the patient's ID.
	Master Patient ID	The unique identifier of the master patient.
	Issuer of Master Patient ID	The unique identifier of the individual that issued the master patient's ID.
	Birth Date	Birth date of the patient.
	Sex	Sex of the patient.
	Is Deceased	Specifies whether the patient is deceased.
	Time of Death	The time of death of the patient.
	Do Not Delete	Specifies whether you want to prevent other users from deleting the patient.

## Edit or delete an event

To edit or delete an event, complete the following steps.

1. In **Worklist**, in the left pane, click **Events**.
2. In the right pane, in the row containing the event that you want to edit or delete, do one of the following:
  - To delete the event, click **Delete** and in the **Delete Event** dialog box, click **Delete Event** to confirm.
  - To edit the event, click **Edit**, edit the event as needed, and then click **Save Changes**. For a description of the available settings, see “Event settings”.

## Event settings

The following table describes the settings that are available when you edit an event. See “Edit or delete an event” for the configuration procedure.

Property Type	Property Value	Property Description
Patient Information	Medical Alerts	Description of patient's medical conditions to which the physician and staff should be alerted.
	Contrast Allergies	Description of patient's reactions or allergies to the contrast agent.
	Patient State	Description of the patient's current state.
	Patient Weight	The weight of the patient.
	Pre Medication	The medication to be administered at the beginning of the procedure.
	Pregnancy Status	Describes the pregnancy state of the patient.
Study Information	Master Accession Number	The master accession number for the procedure.
	Issuer of Master Acc No	The unique identifier of the individual that issued the master accession number for the procedure.
	Study Instance UID	The unique identifier for the study.
Visit	Admissions ID	The unique identifier of the visit as assigned by the healthcare provider.
	Current Patient Location	Description of the current known location of the patient.

<b>Property Type</b>	<b>Property Value</b>	<b>Property Description</b>
	Special Needs	Description of any special needs that the patient may have.
	Transportation Agreement	Mode of transportation of the patient to the location of the procedure.
Requested Procedure	Requested Procedure ID	The unique identifier of the requested procedure.
	Requested Procedure Step ID	The unique identifier of the requested procedure step.
	Requested Procedure Description	The institution-generated administrative description or classification of the requested procedure.
	Requesting Physician	The name of the physician requesting the procedure.
	Referring Physician	The name of the patient's referring physician.
	Requested Contrast Agent	The name of the contrast agent requested for use in the procedure.
	Requested Procedure Priority	The urgency type of the requested procedure.
Scheduled Procedure	Station Name	The name of the modality on which the procedure is scheduled to be performed.
	Station AE Title	The AE title of the modality on which the procedure is scheduled to be performed.
	Scheduled Procedure Start Date	The scheduled start date of the procedure.
	Performing Physician	The name of the physician performing the procedure.
	Accession Number	The unique identifier of the imaging service request for this requested procedure.
	Issuer of Accession Number	The unique identifier of the individual that issued the accession number.

Property Type	Property Value	Property Description
	Procedure Step Description	The institution-generated administrative description or classification of the procedure step.

## Add a worklist connection

To add a worklist connection, complete the following steps.

1. In **Worklist**, in the left pane, click **Configuration**.
2. On the **Configuration** page, click **Add a new worklist configuration**.
3. In the **Connection Details** dialog box, in the **Name** box, type a name for the worklist connection.
4. Select the appropriate worklist connection type, **DICOM Modality Worklist** or **Acuo Semantix**.

**Note:** If you select Acuo Semantix, you only need to configure the settings under Force Patients to REM.

5. In the **Hostname** and **Port** boxes, type the host name and port number of the destination device.
6. In the **Called AE** and **Calling AE** boxes, type the called AE title and calling AE title of the destination device.
7. In the **Query Days Back** and **Query Days Forward** boxes, type the number of days you want the query to search before and after the current date.
8. In the **Request Timer** box, type how often, in seconds, you want the C-FIND to run.
9. In the **Cleanup After** box, type how long, in days, you want query results to be kept. A value of 0 (Zero) will keep results forever.
10. In the **Patient Name** or **Patient ID** boxes, type either the designated fake patient name or ID, such as Doe^Jane^^^ or ER1234.

**Note:** Emergency care patients sometimes arrive into a care facility where the patient's name and ID are unknown. Configuring this step provides a way to force these types of patients into the REM by configuring a designated name or ID to identify them.

11. Click **Save Changes**.

## Edit or delete a worklist connection

To edit or delete a worklist connection, complete the following steps.

1. In **Worklist**, in the left pane, click **Configuration**.
2. In the right pane, in the row containing the worklist that you want to edit or delete, do one of the following:
  - To delete the worklist, click **Delete**.
  - To edit the worklist, click **Edit**, edit the worklist as needed, and then click **Save Changes**.

# XDS Explorer

## What is XDS Explorer?

XDS Explorer is a Acuo Admin Portal module that administrators can use to browse documents in an XDS repository and view associated metadata. Within each repository, you can select which master patient index to use in a search. XDS Explorer also provides a VNA search function to view VNA records alongside XDS document metadata.

## Configure XDS

**Prerequisite** See 'XDS Explorer Requirements' in the technical specifications.

Before you can use XDS Explorer, you must configure an XDS repository in the Administration module. To configure XDS, complete the following steps.

1. In **Administration**, select the **XDS Configuration** tab.
2. In the left pane, click **Add New Acuo Registry**.
3. In the **Add New Acuo Registry** dialog box, type the name, URL, and domain of the registry and click **Save**.
4. In the right pane, click **Add New Acuo MPI**.
5. In the **Add New Acuo MPI** dialog box, type the name and URL of the medical patient index and click **Save**.

## About PDQ and registry searches

### PDQ search

The patient demographic query (PDQ) is the primary means of searching the registry of an XDS repository. Acuo Admin Portal uses the patient demographic criteria that you enter to search the enterprise master patient index (EMPI) for a local ID, and then uses the local ID to find the mapped affinity domain or global ID. Acuo Admin Portal retrieves the matching documents from the registry by using that affinity domain or global ID.

### Registry search

The registry search is less accurate than the PDQ search, but it is very useful for finding erroneous patient data. In a registry search, you enter the same demographic search criteria as in a PDQ search. However, instead of searching the EMPI for an ID, Acuo Admin Portal searches the registry directly for matching criteria in document metadata.

## Perform a PDQ search

To find a patient document in an XDS repository, you can perform a PDQ search. To perform a PDQ search, complete the following steps.

1. In **XDS Explorer**, select a system.
2. On the menu bar, on the **Acuo Registry** tab, select a registry to search.
3. On the menu bar, click **PDQ Search**.
4. In the **PDQ Patient Search** dialog box, complete the following substeps.

1. In the **Acuo MPI** list, select a master patient index.
2. Under **Search Criteria**, type or select one or more criteria and click **Search**.
3. In the **PDQ Results** list, select a patient record and click **Retrieve Documents**.
5. To view document metadata, in the **Registry Documents** list, on the **Actions** menu, select **View**.

## Perform a registry search

You can search for patient demographics in document metadata and return the matching documents. To perform a registry search, complete the following steps.

1. In **XDS Explorer**, select a system.
2. On the menu bar, click **Registry Search**.
3. In the **Registry Search** dialog box, enter search criteria and click **Search**.
4. In the **Registry Results** list, select a patient record and click **Retrieve Documents**.

## Perform an MPQ search

To find all patients in the registry matching non-demographic criteria, you can perform a multiple patient query (MPQ) search. To perform an MPQ search, complete the following steps.

1. In **XDS Explorer**, select a system.
2. On the menu bar, click **Registry Search**.
3. In the **Registry Search** dialog box, select the **Multiple Patient Search (MPQ)** check box.
4. Enter search criteria and click **Retrieve Documents**.

## Perform a VNA search

Although the main purpose of XDS Explorer is to search repositories, you can also use it to search the VNA. To perform a VNA search, complete the following steps.

**Note:** The VNA search option is only available if you do not have an ACUO EMPI configured.

1. In **XDS Explorer**, select a system.
2. On the menu bar, click **VNA Search**.
3. In the **VNA Patient Search** dialog box, complete the following substeps.
  1. In the **Sources** list, select a database.
  2. Under **Search Criteria**, type or select one or more criteria and click **Search**.
  3. In the **VNA Results** list, select a patient record and click **Retrieve Documents**.

# System Management

## What is System Management?

System Management is a Acuo Admin Portal module that administrators can use to view and management the health of the Acuo Admin Portal system. The module provides three areas of configuration: application servers, services, and performance.

### Application Servers

Displays a list of application servers, and enables you to start, stop, and restart services associated to those servers.

### Services

Displays a list of Acuo services and enables you restart all those services collectively.

### Performance

Displays a list of application servers and databases, and enables you to view system performance.

## Start a service

To start a service, complete the following steps.

1. In **System Management**, in the left pane, click **Application Servers**.
2. On the **Application Servers** page, locate the service you want to start and click **Manage Service**.
3. In the **Manage Service** dialog box, click **Start**. The system starts the service.

## Stop a service

To stop a service, complete the following steps.

1. In **System Management**, in the left pane, click **Application Servers**.
2. On the **Application Servers** page, locate the service you want to stop and click **Manage Service**.
3. In the **Manage Service** dialog box, click **Stop**. The system stops the service.

## Restart a service

To restart a service, complete the following steps.

1. In **System Management**, in the left pane, click **Application Servers**.
2. On the **Application Servers** page, locate the service you want to restart and click **Manage Service**.
3. In the **Manage Service** dialog box, click **Restart**. The system restarts the service.

## Configure Application Server with a load balancer

To configure the application server with a load-balancer, complete the following steps.

1. In **System Management**, in the left pane, click **Application Servers**.
2. On the **Application Servers** page, locate the application server you want to add to or remove from the load balancer.
3. Beside **Is In Load Balancer**, click its current state, **Yes** or **No**.
4. In the **Application Server** dialog box, select **Is In Load Balancer** to add the application server to the load balancer or deselect **Is In Load Balancer** to remove the application server from the load balancer.
5. Click **Save Changes**.

## View a list of services

To view a list of services, complete the following steps.

1. In **System Management**, in the left pane, click **Services**.
2. On the **Services** page, the system lists all the services currently available for each application server instance.
3. Optional. Click **Refresh** to update the information.

## Restart all services

To restart all services, complete the following steps.

1. In **System Management**, in the left pane, click **Services**.
2. On the **Services** page, click **Restart All Services**.
3. In the **Restart All Services** dialog box, click **Restart All Services** to confirm.
4. In the **Restart All Services** dialog box, click **Close**. The system restarts all services.

## View system performance

You can monitor system performance by displaying a summary of application servers and databases. To display this list, complete the following steps.

1. In **System Management**, in the left pane, click **Performance**.
2. On the **Performance** page, the system lists a summary of performance information for application servers, processes, disk I/O, SQL servers and SQL databases.
3. Optional. Click **Refresh** to update the information.

# Lifecycle Management

## What is lifecycle management?

Lifecycle Management is a Acuo Admin Portal module that administrators can use to help manage data volume. Lifecycle management refers to the scheduled moving or deleting of records from a system according to retention policies. Lifecycle management requires two tasks, creating a retention policy and scheduling a task that executes the retention policy. The Lifecycle Management module includes a sub-module for each of these two tasks.

## Create a retention policy

You can create a retention policy to specify when Acuo Admin Portal should delete specific records. Acuo Admin Portal does not execute the policy unless you schedule it in Task Scheduling. To create a retention policy, complete the following steps.

1. In **Retention Policies**, click **Create New Policy**.
2. In the **Create New Policy** dialog box, type a name for the policy and enter one or more of the following criteria.
  - Under **Patient Birthdate**, **Study Date**, or **Insert Date**, select the check box and specify a number of days to retain the record after the birth date of the patient, the study date, or the insert date.
  - Under **Institution Name**, select the check box and type the name of an institution. The policy only applies to records from this institution.



- In the **Modalities** list, select one or more modalities, and then click anywhere outside the list. The policy only applies to records from the modalities you select.
- 3. In the **Action** list, select one of the following actions to perform when the policy is executed.
  - Select **Delete Images** if you want the system to delete all images.
  - Select **Move Images to Recycle Bin** if you want the system to move all images to the recycle bin.
  - Select **Compress and Move Images to an ILM Database** if you want the system to compress all images and move them to an ILM database.
- 4. Optional. To limit the scope of the policy by destination, complete the following substeps.
  1. Clear the **Apply to all destinations** check box.
  2. In the **Destinations** list, select one or more destinations and click anywhere outside the list.
- 5. Click **Create New Policy**.

## Schedule execution of a retention policy

After you create a retention policy you must schedule a task to execute it. Tasks apply to all policies of a specified destination. The task that you create runs in the background as a Windows scheduled task. To schedule execution of a retention policy, complete the following steps.

1. In **Task Scheduling**, click **Create New Scheduled Task**.
2. In the **Create New Task** dialog box, complete the following substeps.
  1. To assign policies to the task, in the **Acuo Destination** list, select a destination.
  2. In the **Frequency** list, select an execution frequency of **Daily**, **Weekly**, or **Monthly**.
  3. Under **Run Time**, specify the time of day to execute the task.
  4. In the **Task Username** and **Task Password** boxes, type credentials for the domain or machine on which the task runs.
3. Click **Create New Task**.

## Auditing

### About Auditing

Acuo Admin Portal maintains four types of audit logs that can be viewed in the following modules.

#### Acuo Audits

Shows access to patient data to help you track who viewed which patients.

#### IHE Audits

Shows event logs conforming to IHE standards.

#### User Audits

Shows detailed activity by user.

#### Tasks

Shows audit information sent to external systems only if AcuoMed is configured to do so.

## View an audit log

You can view audit logs of server activity. To view an audit log, complete the following steps.

1. On the **Auditing** module home screen, select **Acuo Audits, IHE Audits, User Audits, or Tasks**.
2. Optional. To customize the log, see *Customize a list*.

## System Configuration

### What is System Configuration?

System Configuration is a Acuo Admin Portal module that administrators can use to edit system configuration. The module provides three areas of configuration: batch processing, IPID domains, and SOPs.

#### Batch Processing

Displays a list of application servers and the associated batch processing options, and enables you to edit the processing options for each application server.

#### IPID Domains

Displays a list of IPID domains, and enables you to view, delete, and edit those domains.

#### SOPs

Displays a list of SOPs, and enables you to edit or delete those SOPs.

### Edit batch processing options

To edit batch processing options, complete the following steps.

1. In **System Configuration**, in the left pane, click **Batch Processing Options**.
2. On the **Batch Processing Options** page, locate the application server you want to update and click **Edit Processing Options**.
3. In the **Edit Processing Options** dialog box, edit the processing options as needed, and then click **Save Changes**. For a description of the available options, see “Batch processing options”.

### Batch processing options

The following table describes the settings that are available when you edit a batch process. See “Editing batch processing options” for the configuration procedure.

Property Value	Property Description
Maximum Job	The number of batch store jobs that can be running at one time.
Keep Completed Jobs	The number of days that jobs will continue to appear in the Batch Store Manager job summary.

Property Value	Property Description
Associated Timeout	The amount of time, in seconds, an association with the DICOM device will stay open after the transmission of an image is complete.
Default Job Priority	The default priority level of the batch store job.
Idle Time Before Group Completes	The amount of time, in seconds, the queue is idle before the next group starts.

## Create an IPID domain

To create an IPID domain, complete the following steps.

1. In **System Configuration**, in the left pane, click **IPID Domains**.
2. On the **IPID Domains** page, click **Create New IPID Domain**.
3. In the **Create New IPID Domain** dialog box, complete the following sub-steps:
  1. In the **Name** field, type a name for the IPID domain.
  2. In the **Issuer of Patient ID** field, type the unique identifier of the assigning authority that issued the patient ID.
  3. In the **Assigning Authority** field, type the name of the Assigning Authority that issued the patient ID.
4. Click **Save Changes**.

## Edit or delete an IPID domain

To edit or delete an IPID domain, complete the following steps.

1. In **System Configuration**, in the left pane, click **IPID Domains**.
2. In the right pane, in the row containing the IPID domain that you want to edit or delete, do one of the following:
  - To delete the IPID domain, click **Delete** and in the **Delete IPID Domain** dialog box, click **Delete IPID Domain** to confirm.
  - To edit the IPID domain, click **Edit**, edit the IPID domain as needed, and then click **Save Changes**.

## Create an SOP

To create an SOP, complete the following steps.

1. In **System Configuration**, in the left pane, click **SOPs**.
2. On the **SOPs** page, click **Create new SOP**.
3. In the **Create new SOP** dialog box, complete the following sub-steps:

1. In the **SOP UID** field, type the unique identifier for the SOP.
2. In the **SOP UID Description** field, type a description for the SOP.
3. From the **SOP Group** list, select an SOP group for this SOP.
4. Click **Save Changes**.

## Edit or delete an SOP

To edit or delete an SOP, complete the following steps.

1. In **System Configuration**, in the left pane, click **SOPs**.
2. In the right pane, in the row containing the IPID domain that you want to edit or delete, do one of the following:
  - To delete the SOP, click **Delete** and in the **Delete SOP** dialog box, click **Delete SOP** to confirm.
  - To edit the SOP, click **Edit**, edit the SOP name, and then click **Save Changes**.