



Gear View Basic
DICOM Conformance Statement
ENG-GVB-DCS-REVC

Date: July 2025

Documentation Notice

Information in this document is subject to change without notice. The software described in this document is furnished only under a separate license agreement and may only be used or copied according to the terms of such agreement. It is against the law to copy the software except as specifically allowed in the license agreement. This document or accompanying materials may contain certain information which is confidential information of Hyland Software, Inc. and its affiliates, and which may be subject to the confidentiality provisions agreed to by you.

Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright law, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the express written permission of Hyland Software, Inc. or one of its affiliates.

Hyland, HXP, OnBase, Alfresco, Nuxeo, and product names are registered and/or unregistered trademarks of Hyland Software, Inc. and its affiliates in the United States and other countries. All other trademarks, service marks, trade names and products of other companies are the property of their respective owners.

© 2025 Hyland Software, Inc. and its affiliates.

The information in this document may contain technology as defined by the Export Administration Regulations (EAR) and could be subject to the Export Control Laws of the U.S. Government including for the EAR and trade and economic sanctions maintained by the Office of Foreign Assets Control as well as the export controls laws of your entity's local jurisdiction. Transfer of such technology by any means to a foreign person, whether in the United States or abroad, could require export licensing or other approval from the U.S. Government and the export authority of your entity's jurisdiction. You are responsible for ensuring that you have any required approvals prior to export.

Revision History

Date	Revision	Authors	Description
17-NOV-2010	A	Chris Barnett	Created.
01-NOV-2012	B	Marilyn Palowitch	Corrected part number, removing use of product version number in the part number; minor document formatting improvements.
02-OCT-2013	C	Chris Barnett	Updated for version 2.1 Added Breast Tomo
10-JUN-2016	A	Cheryl Hawkins	Rebrand Manufacturer name from Pacsgear Inc to Lexmark Enterprise Software, LLC. Update and logo. Updated copyright and new Product ID. Changes to Product Name from GEARView to Gear View – excluding current AE name of GEARView.
09-SEP-2023	B	Tammy Matthews	Rebranded for Hyland. Minor formatting changes.
09-APR-2025	C	Tammy Matthews	Rebranded for Hyland 2025 initiative.

Table of Contents

Documentation Notice.....	2
Revision History	3
List of Figures.....	5
Introduction	6
Implementation Model.....	6
Application Data Flow Diagram	6
Functional Definitions of Application Entity (AE)	6
Sequencing of Real-World Activities.....	6
AE Specifications	6
Gear View Basic AE Specifications	6
Functional Definition of Application Entity	Error! Bookmark not defined.
<i>Sequencing of Real-World Activities – Read Images from Media.....</i>	6
Communication Profiles.....	8
Supported Communication Stacks.....	8
TCP/IP Stack	8
Physical Media Support.....	8
Extension/Specialization/Privatization.....	9
Configuration	9
Extended Character Sets	9

List of Figures

Figure 1.	Gear View Basic Implementation Model.....	6
Figure 2.	Gear View Basic Application Profiles.....	6
Figure 3.	Gear View Basic Supported SOP Classes	8

Introduction

This conformance statement is designed to communicate technical information regarding the Gear View Basic product and its compliance to the DICOM 3.0 standard. Gear View Basic is a DICOM viewer that lets physicians and patients view medical images and related results.

Implementation Model

Application Data Flow Diagram




Figure 1. Gear View Basic Implementation Model

Functional Definitions of Application Entity (AE)

Not applicable.

Sequencing of Real-World Activities

Not applicable.

AE Specifications

Functional Definition of AE

Sequencing of Real-World Activities – Read Images from Media

The Gear View Basic AE provides standard conformance to the following DICOM 3.0 Interchange option for Media Storage service class with the following profiles and roles.

Application Profiles	Real World Activity	Role	SC Option
STD-GEN-CD	Read a CD	FSR	Interchange
STD-GEN-DVD	Read a DVD	FSR	Interchange

Figure 2. Gear View Basic Application Profiles

Supported SOP Classes

The Gear View Basic AE reads the following SOP classes from media.

SOP Class Name	SOP Class UID
Breast Tomosynthesis Image Storage	1.2.840.10008.5.1.4.1.1.13.1.3
Computed Radiography Image Storage	1.2.840.10008.5.1.4.1.1.1
Digital X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.1
Digital X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.1.1
Digital Mammography Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.2
Digital Mammography Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.2.1
Digital Intra Oral X-Ray Image Storage – For Presentation	1.2.840.10008.5.1.4.1.1.1.3
Digital Intra Oral X-Ray Image Storage – For Processing	1.2.840.10008.5.1.4.1.1.1.3.1
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1
X-Ray Angiographic Image Storage	1.2.840.10008.5.1.4.1.1.12.1
Enhanced XA Image Storage	1.2.840.10008.5.1.4.1.1.12.1.1
X-Ray Radiofluoroscopic Image Storage	1.2.840.10008.5.1.4.1.1.12.2
Positron Emission Tomography Image Storage	1.2.840.10008.5.1.4.1.1.128
CT Image Storage	1.2.840.10008.5.1.4.1.1.2
Enhanced CT Image Storage	1.2.840.10008.5.1.4.1.1.2.1
Nuclear Medicine Image Storage	1.2.840.10008.5.1.4.1.1.20
US Multi-Frame Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.3
US Multi-Frame Image Storage	1.2.840.10008.5.1.4.1.1.3.1
MR Image Storage	1.2.840.10008.5.1.4.1.1.4
Enhanced MR Image Storage	1.2.840.10008.5.1.4.1.1.4.1
RT Image Storage	1.2.840.10008.5.1.4.1.1.481.1
NM Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.5
US Image Storage (Retired)	1.2.840.10008.5.1.4.1.1.6
Ultrasound Image Storage	1.2.840.10008.5.1.4.1.1.6.1

Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7
Multi-frame Grayscale Byte Secondary Capture Image Storage	1.2.840.10008.5.1.4.1.1.7.2
Multi-frame True Color Secondary Capture Image	1.2.840.10008.5.1.4.1.1.7.4
VL Endoscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.1
VL Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.2
VL Slide Coordinates Microscopic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.3
VL Photographic Image Storage	1.2.840.10008.5.1.4.1.1.77.1.4
Ophthalmic Photography 8-Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.1
Ophthalmic Photography 16-Bit Image Storage	1.2.840.10008.5.1.4.1.1.77.1.5.2
Basic Text SR†	1.2.840.10008.5.1.4.1.1.88.11
Enhanced SR Storage†	1.2.840.10008.5.1.4.1.1.88.22
Comprehensive SR Storage†	1.2.840.10008.5.1.4.1.1.88.33
Grayscale Softcopy Presentation State Storage SOP Class*	1.2.840.10008.5.1.4.1.1.11.1
Mammography CAD SR Transfer†	1.2.840.10008.5.1.4.1.1.88.50
Key Object Selection Document†	1.2.840.10008.5.1.4.1.1.88.59

Figure 3. Gear View Basic Supported SOP Classes

* - This modality will display a "modality not currently supported" message.

† - This modality will display as plain text

Communication Profiles

Supported Communication Stacks

Not applicable.

TCP/IP Stack

Not applicable.

Physical Media Support

Not applicable.

Extension/Specialization/Privatization

Not applicable.

Configuration

Not applicable.

Extended Character Sets

Not applicable.