



Media Writer  
Alternate Viewer Installation Guide

DOC-MW-AV-IN-EN-REVB

# Documentation Notice

Information in this document is subject to change without notice. The software described in this document is furnished only under a separate license agreement and may only be used or copied according to the terms of such agreement. It is against the law to copy the software except as specifically allowed in the license agreement. This document or accompanying materials may contain certain information which is confidential information of Hyland Software, Inc. and its affiliates, and which may be subject to the confidentiality provisions agreed to by you.

Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright law, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the express written permission of Hyland Software, Inc. or one of its affiliates.

Hyland, HXP, OnBase, Alfresco, Nuxeo, and product names are registered and/or unregistered trademarks of Hyland Software, Inc. and its affiliates in the United States and other countries. All other trademarks, service marks, trade names and products of other companies are the property of their respective owners.

© 2023 Hyland Software, Inc. and its affiliates.

The information in this document may contain technology as defined by the Export Administration Regulations (EAR) and could be subject to the Export Control Laws of the U.S. Government including for the EAR and trade and economic sanctions maintained by the Office of Foreign Assets Control as well as the export controls laws of your entity's local jurisdiction. Transfer of such technology by any means to a foreign person, whether in the United States or abroad, could require export licensing or other approval from the U.S. Government and the export authority of your entity's jurisdiction. You are responsible for ensuring that you have any required approvals prior to export.



## Regulatory and compliance

**Australian Sponsor:**

Emergo Australia  
Level 20 Tower II  
Darling Park 201 Sussex Street  
Sydney, NSW 2000 Australia



**Manufacturer**  
**Hyland Software, Inc.**  
12919 Earhart Ave  
Auburn, CA 95602

**Support**

Telephone: 1-844-535-1404  
Web: <https://community.hyland.com>

## General warnings and precautions



To ensure patient safety and adhere to HIPAA Privacy Rule involving patient health information, the Healthcare provider shall assign an unique Accession Number and Medical Record Number (MRN) to identify individual patient's health record affiliated with the patient to prevent erroneous medical data and duplicate identification with other patients.

## Symbols glossary



3082 - Manufacturer: Indicates the manufacturer.



2497 - Date of Manufacture: Indicates the date when the medical device was manufactured.



2493 - Catalog number: Indicates the manufacturer's catalog number so the medical device can be identified.



2498 - Serial number: Indicates the manufacturer's serial number so that a specific medical device can be identified.



1641 - Consult instructions for use: Indicates the need for the user to consult the instructions for use.



0434A - Caution: Indicates the need for the user to consult the instructions for use for important information such as warnings and cautions.



2610 - Patient number: Indicates a unique number associated with an individual patient.



Quantity: Indicates the number of units.



Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

# Table of Contents

Documentation Notice ..... 2

Regulatory and compliance ..... 4

General warnings and precautions ..... 4

Symbols glossary ..... 5

About viewers for Media Writer ..... 7

Install and alternate viewer ..... 7

Edit the viewer configuration file ..... 7

## About viewers for Media Writer

Media Writer comes with a viewer called Gear View Basic that you can include on burned media so that patients can view DICOM studies. If you prefer to include a different viewer, you must first install it on the Media Writer server computer and then add it as part of your Media Writer installation. The alternate viewers that are available are eFilm 3.4, MIMVista, and ThreePalm.

## Install and alternate viewer

To install an alternate viewer, complete the following steps.

### Prerequisite

Obtain the USB flash drive containing the viewer files and configuration files from your Lexmark sales representative.

1. On the server computer for Media Writer, create a new folder called *Original* under the *C:/Program Files (x86)/Pacsgear/Media Writer/Service/Config* folder.
2. Move the default *viewers.config* file from the *Config* folder into the *Original* folder.
3. Copy the *xyz\_viewers.config* file from the USB flash drive to *Media Writer/Service/Config*.

**Note:** xyz represents efilm, mim, or tsp depending on which viewer you want to install.

4. Rename the *xyz\_viewers.config* file to *viewers.config*.
5. Copy the appropriate Viewer folder from the USB flash drive to *../Media Writer/Service/Viewers*.
6. If installing the ThreePalm Mammoviewer, copy *postinstall.exe* from the USB flash drive to *../Media Writer/Service*.
7. Close all windows, eject the USB flash drive, and then restart the computer.

## Edit the viewer configuration file

After installing the alternate viewer and configuration file, you must edit the configuration file of the new viewer to make it available from within Media Writer. To edit the viewer configuration file, complete the following steps.

1. In **Control Panel**, click **Administrative Tools**, and then double-click **Services**.
2. Select the Media Writer service in the list, and then click **Stop**.
3. Open the *../Media Writer/Service/Config/viewers.config* file in Notepad.
4. At the bottom of the file, add the text below that corresponds to the viewer you installed.

Viewer	Add to configuration file
eFilm	<pre>&lt;VIEWER&gt; &lt;NAME-&gt;eFilm&lt;/NAME&gt; &lt;PATH&gt;viewer/eFilm&lt;/PATH&gt; &lt;VIEWER_EXE_FILE&gt;efilmLt.exe&lt;/VIEWER_EXE_FILE&gt; &lt;/VIEWER&gt;</pre>

Viewer	Add to configuration file
3Palm Mammo	<pre> &lt;VIEWER&gt; &lt;NAME-&gt;MAMMO TPS&lt;/NAME&gt; &lt;PATH&gt;<i>[relative path to the viewer folder]</i>&lt;/PATH&gt; &lt;VIEWER_EXE_FILE&gt;<i>[full name of the viewer executable file]</i>&lt;/VIEWER_EXE_ FILE&gt; &lt;/VIEWER&gt; &lt;VIEWER_LICGEN_EXE&gt;postinstall.exe&lt;/VIEWER_LICGEN_EXE &gt; </pre>
MIMVista	<pre> &lt;VIEWER&gt; &lt;NAME-&gt;MIMVISTA&lt;/NAME&gt; &lt;PATH&gt;<i>[relative path to the viewer folder]</i>&lt;/PATH&gt; &lt;VIEWER_EXE_FILE&gt;<i>[full name of the viewer executable file]</i>&lt;/VIEWER_EXE_ FILE&gt; &lt;/VIEWER&gt; &lt;VIEWER_LICGEN_EXE&gt;mvprep.exe&lt;/VIEWER_LICGEN_EXE &gt; </pre>

5. Click **Save** and close Notepad.
6. Restart the Media Writer service