

# PACS Scan Mobile

## Release Notes

Version: 2.0.1

Written by: Product Knowledge, R&D

Date: September 2016

LX-DOC-PSM2.0.1-RN-EN-REVB

**Australian Sponsor**

Emergo Australia  
Level 20 Tower II  
Darling Park 201 Sussex Street  
Sydney, NSW 2000 Australia



**Lexmark Enterprise Software, LLC**

4309 Hacienda Drive, Suite 500  
Pleasanton, CA 94588 USA



**EMERGO EUROPE**  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



© 2016 Lexmark. All rights reserved.

**Lexmark is a trademark of Lexmark International Inc., registered in the U.S. and/or other countries. All other trademarks are the property of their respective owners. No part of this publication may be reproduced, stored, or transmitted in any form without the prior written permission of Lexmark.**

## Table of Contents

<b>PACS Scan Mobile 2.0.x (March 2016)</b> .....	<b>4</b>
Enhancements.....	4
Updates .....	4
<b>PACS Scan Mobile 1.1.x (October 2015)</b> .....	<b>4</b>
Enhancements.....	4
<b>PACS Scan Mobile 1.0.x (April 2015)</b> .....	<b>4</b>
About PACS Scan Server.....	4
About PACS Scan Mobile Administration Console .....	4
Features.....	5

## PACS Scan Mobile 2.0.x (March 2016)

The following enhancements and updates are included in this release.

### Enhancements

- You can now provision a mobile device running PACS Scan Mobile by using a QR code.
- You can now log in with Touch ID.
- PACS Scan Mobile can now send studies to XDS repositories.
- Users can now launch PACS Scan Mobile with encrypted arguments from other iOS apps such as Epic Haiku and Canto.

### Updates

- The application now checks the version and rebuilds the database if necessary.
- Changed the DICOM command on the Settings menu to Connections.

## PACS Scan Mobile 1.1.x (October 2015)

The following enhancements and updates are included in this release.

### Enhancements

- PACS Scan Mobile is now available on iTunes.
- Administrators can now provision devices for users by setting the IP address and SSL options.
- You can now include a backslash in an MRN for patient entry and lookup.

## PACS Scan Mobile 1.0.x (April 2015)

PACS Scan Mobile is an app for mobile devices that physicians and authorized staff can use to capture images and videos easily. You can use the built-in camera of the mobile device to capture images or videos of your patient and then attach them to patient studies. When you are finished capturing images and videos, you must send them to a selected query source for storage, but you have the option of editing the images first.

### About PACS Scan Server

PACS Scan Server is a desktop program that controls the flow of data between DICOM or related systems such as PACS, archives, Worklist servers, modalities, EHR/EMR, and RIS, and any clients that connect to it such as the PACS Scan Mobile app.

### About PACS Scan Mobile Administration Console

PACS Scan Mobile Administration Console (PSM Console) is a browser-based interface to PACS Scan Server that administrators can use to configure server settings, and settings for all of the client mobile devices that connect to the server. These include settings for licensing, DICOM querying and sending, users, and logging.

## Features

PACS Scan Mobile offers the following features.

- Query and display patient demographics from multiple sources, or create new patient data.
- Capture clinical images and videos, attach to patient studies, and send to DICOM archives.
- Edit demographics and media.
- Manage mobile device configuration and user access from a web console.
- Create department-specific configurations.
- Make both DICOM and Worklist queries, and send HL7® messages.
- View DICOM Trace data, server logs, and audit logs.