



# PACS Scan Mobile

## Release Notes


DOC-PSM2.2.1-RN-EN-REVA

## Regulations and Compliance

**Australian Sponsor:**

Emergo Australia  
Level 20 Tower II  
Darling Park 201 Sussex Street  
Sydney, NSW 2000 Australia

**EU Sponsor:**



**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



**Hyland LLC**

4309 Hacienda Drive, Suite 500  
Pleasanton, CA 94588 USA

**Tel: 1-844-535-1404**

**Email: [TS\\_PACSGEAR@hyland.com](mailto:TS_PACSGEAR@hyland.com)**

**© 2018 Hyland. Hyland and the Hyland logo are trademarks of Hyland LLC, registered in the United States and/or other countries. All other trademarks are the property of their respective owners.**

# Table of Contents

<b>Regulations and Compliance</b> .....	<b>2</b>
<b>What is PACS Scan Mobile?</b> .....	<b>4</b>
<b>PACS Scan Mobile 2.2.1 (December 2018)</b> .....	<b>4</b>
Features .....	4
<b>PACS Scan Mobile 2.2.0 (March 2018)</b> .....	<b>4</b>
Features .....	4
Updates .....	4
<b>PACS Scan Mobile 2.1.1 (November 2017)</b> .....	<b>4</b>
Features .....	4
<b>PACS Scan Mobile 2.1.0 (September 2017)</b> .....	<b>5</b>
Features .....	5
Product enhancements .....	5
<b>PACS Scan Mobile 1.1.x (October 2015)</b> .....	<b>5</b>
Product enhancements .....	5
<b>PACS Scan Mobile 1.0.x (April 2015)</b> .....	<b>5</b>
Features .....	5

## What is PACS Scan Mobile?

PACS Scan Mobile is an app for iOS and Android mobile devices that physicians and authorized staff can use to capture images and videos. You can use the built-in camera of the mobile device to capture images or videos of your patient and then attach them to patient studies. When you are finished capturing images and videos you can edit them, send them to the selected source for storage, and optionally send an HL7 message to an HL7 device. You cannot save images and videos on your mobile device.

PACS Scan Mobile app is a client of Pacsgear Core Server, which processes DICOM and non-DICOM medical data, and handles user accounts, security, client configuration, and the flow of data between modalities and destination DICOM systems.

## PACS Scan Mobile 2.2.1 (December 2018)

This release includes the following features and updates. The features require PACSgear Core Server 4.1 or later.

### Features

- PACS Scan Mobile is now available for Android phones and tablets. The minimum operating system requirement is Android 6.0 Marshmallow.

#### Note:

Currently, Android PACS Scan Mobile can only import one image at a time.

## PACS Scan Mobile 2.2.0 (March 2018)

This release includes the following features and updates. The features require PACSgear Core Server 3.6.0.

### Features

- You can now annotate images with a fingertip or stylus.
- You can now retrieve patient records by scanning a barcode.

### Updates

- PACS Scan Mobile now correctly sets the MIME type when sending video to XDS destinations.
- Resolved an issue with body mapping when sending to XDS destinations.

## PACS Scan Mobile 2.1.1 (November 2017)

This release includes the following features.

### Features

- When applications launch PACS Scan Mobile with a URL, they can now pass the user name in that URL which is then stored in the Operator DICOM tag (0008,1070) and the AuthorName field in an XDS submission.

## PACS Scan Mobile 2.1.0 (September 2017)

This release includes the following features and product enhancements.

### Features

- You can now append a body part description to studies.
- You can now import photos and videos from the camera roll of the mobile device.
- You can now zoom in and out prior to capturing images.

### Product enhancements

- When sending studies, PACS Scan Mobile now includes the referring physician and additional patient history in the metadata.

## PACS Scan Mobile 1.1.x (October 2015)

This release includes the following product enhancements.

### Product enhancements

- PACS Scan Mobile is now available on iTunes.
- Administrators can now provision devices for users by setting the IP address and SSL options.
- You can now include a backslash in an MRN for patient entry and lookup.

## PACS Scan Mobile 1.0.x (April 2015)

This release includes the following features.

### Features

- Query and display patient demographics from multiple sources, or create new patient data.
- Capture clinical images and videos, attach to patient studies, and send to DICOM archives.
- Edit demographics and media.
- Manage mobile device configuration and user access from a web console.
- Create department-specific configurations.
- Make both DICOM and Worklist queries, and send HL7® messages.
- View DICOM Trace data, server logs, and audit logs.