



PACSGear Enterprise Server Installation Guide

DOC-PSM3.2-IN-EN-REVA

Copyright


Information in this document is subject to change without notice. The software described in this document is furnished only under a separate license agreement and may be used or copied only according to the terms of such agreement. It is against the law to copy the software except as specifically allowed in the license agreement. This document or accompanying materials contains certain information which is confidential information of Hyland Software, Inc. and its affiliates, and which is subject to the confidentiality provisions agreed to by you.

All data, names, and formats used in this document's examples are fictitious unless noted otherwise. Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright law, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the express written permission of Hyland Software, Inc. or one of its affiliates.

Hyland®, Hyland Software®, Hyland Healthcare, and Hyland product names, are registered and/or unregistered trademarks of Hyland Software, Inc. and its affiliates in the United States and other countries. All other trademarks, service marks, trade names and products of other companies are the property of their respective owners.

© 2022 Hyland Software, Inc. and its affiliates. All rights reserved.

Regulatory and compliance

Australian Sponsor: Emergo Australia Level 20 Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia	PACS Scan Mobile Anvisa n°: 80117589010 Responsável Técnico: Luiz Levy Cruz Martins CRF-SP: 42415 IMPORTADOR: Emergo Brazil Import Importação e Distribuição de Produtos Médicos Hospitalares Ltda. Avenida Francisco Matarazzo, 1.752, Salas 502/503, Água Branca, São Paulo-SP, CEP - 05001-200 CNPJ: 04.967.408/0001-98 e-mail: brazilvigilance@ul.com
EU Representative:  EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands	



Manufacturer

Hyland LLC

4309 Hacienda Drive, Suite 500
Pleasanton, CA 94588 USA

Support

Telephone: 1-844-535-1404

Web: <https://community.hyland.com>

General warning and precautions



To ensure patient safety and adhere to HIPAA Privacy Rule involving patient health information, the Healthcare provider shall assign an unique Accession Number and Medical Record Number (MRN) to identify individual patient's health record affiliated with the patient to prevent erroneous medical data and duplicate identification with other patients.

Symbols glossary



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

3082 - Manufacturer: Indicates the medical device manufacturer.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

Authorized representative in European Community: Indicates the authorized representative in the European union. Symbol is accompanied by the name and address of the authorized representative adjacent to the symbol.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2497 - Date of Manufacture: Indicates the date when the medical device was manufactured.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2493 - Catalog number: Indicates the manufacturer's catalog number so the medical device can be identified.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2498 - Serial number: Indicates the manufacturer's serial number so that a specific medical device can be identified.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

1641 - Consult instructions for use: Indicates the need for the user to consult the instructions for use.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

0434A - Caution: Indicates the need for the user to consult the instructions for use for important information such as warnings and cautions.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2610 - Patient number: Indicates a unique number associated with an individual patient.



Medical Devices Directive 2017/745/EC

CE Marking: Indicates that a product complies with applicable European Union regulations.



Quantity: Indicates the number of units.



Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).



ISO/DIS ISO 15223-1 - Medical devices - Symbols to be applied to be used with medical device labels, and information to be supplied.

Indicates the item is a medical device.

Table of Contents

- Copyright 2**
- Regulatory and compliance 3**
- General warning and precautions 3**
- Symbols glossary 4**
- What is PACS Scan Mobile? 7**
- Install PACSgear Enterprise Server 7**
- License PACS Scan Mobile 7**
- Install the PACS Scan Mobile app and provision your device 7**
- Confirm the server address 8**
- Log in to PACS Scan Mobile for the first time 8**
- Log in to PACS Scan Mobile with Touch ID or Face ID 8**

What is PACS Scan Mobile?

The PACS Scan Mobile client is an Android and iOS application for mobile devices that allows you to easily capture images and videos and then send them to patient studies on DICOM archives. With proper authorization, you can query and retrieve a patient study from a worklist or create a new patient record. You can take photos and videos. You can manipulate photos, for example, crop or rotate, and also annotate them prior to sending to other locations. You cannot alter videos prior to sending. You can delete unwanted images and videos prior to sending.

Install PACSgear Enterprise Server

PACS Scan Mobile is included with PACSgear Enterprise Server. Install PACSgear Enterprise Server by following the instructions in the PACSgear Enterprise Server Installation and Setup Guide.

License PACS Scan Mobile



Prerequisite Obtain a license key by sending an email to key-pacsgear@hyland.com. In the body of the email, type the following information: Hospital Name, City, State, Country/Region, Product Name, System ID, and Department/Workstation.

If PACS Scan Mobile was included in your purchase of Pacsgear Enterprise Server when you first installed Pacsgear Enterprise Server, you do not need to license PACS Scan Mobile and you can skip this procedure. If you purchase a new PACS Scan Mobile license for an existing Pacsgear Enterprise Server installation, complete the following steps to update your license.

1. From a web browser, go to <http://localhost/pgapps>.

Note: If not logging in from the server computer, substitute localhost with the IP address of the server computer.

2. In PACSgear Enterprise Server, on the **Help** menu, select **About**.
3. In the list of products, next to PACS Scan Mobile, click **License key**.
4. Click **Change** and enter your license key.
5. Click **Save**.

Result After licensing PACS Scan Mobile, the Login page opens.

Install the PACS Scan Mobile app and provision your device

Prerequisite Install and license Pacsgear Enterprise Server.

To install the PACS Scan Mobile app and provision your mobile device, complete the following steps.

1. In the iTunes store, install and open the PACS Scan Mobile app.
2. On a PC, log in to Pacsgear Enterprise Server.
3. On the PACS Scan Mobile app, tap the QR code icon.
The camera of the mobile device turns on.
4. Point the camera at the QR code on the PC screen until the **Settings Updated** message appears.

Confirm the server address

To confirm the server address, complete the following steps.

Your PACS Scan Mobile app registers the address of its Pacsgear Enterprise Server instance. If you have questions about the server address, contact your support representative.

1. On the mobile device, launch the app.
2. To view the server address, on the login screen, tap the information icon.
3. Tap anywhere to return to the login screen.

Log in to PACS Scan Mobile for the first time

To log in to PACS Scan Mobile for the first time, complete the following steps.

1. On the log in screen, type your user name and password and tap **Log In**.
2. Optional. To save your thumbprint, click **Yes**, and then press your thumb on the thumbprint reader.

Note: This step is only available if you registered your thumbprint when you set up your mobile device.

3. Optional. If your device supports **Face ID**, click **Yes** to enable the **Face ID** log in.

Log in to PACS Scan Mobile with Touch ID or Face ID

Prerequisite Your mobile device must include a Touch ID or Face ID function.

If you saved your thumbprint the first time you logged in, you can log in again using Touch ID. To log in using Touch ID, complete the following steps.

1. On the log in screen, tap the thumbprint icon.
2. Press your thumb on the thumbprint reader.
3. If you enabled **Face ID** the first time you logged in, you can log in again using your **Face ID**. To log in using your face, complete the following sub-steps.
 1. On the log in screen, tap the **Face** icon.
 2. Perform a **Face ID** unlock.