



Media Writer

IHE Integration Statement

LX-ENG-MW-IHEIS-EN-REVA

Revision history

Author	Date	Comment
Chris Barnett	12/07/2007	Initial Version, 1.0 REVA
Chris Barnett	7/17/2008	Updated for MW 2.0 REVA
Todd Wysuph	10/12/2009	Minor editing and updating, 2.0 REVB
Todd Wysuph	05/15/2018	Rebranding and updating for MW 5.1.9, REVA

Introduction

This IHE integration statement describes the conformance of this PACSgear product with the IHE Technical Framework. The IHE Technical Framework defines a subset of the functional components of the healthcare enterprise called *IHE Actors*, and specifies their interactions in terms of a set of coordinated standards-based transactions. These transactions are organized into functional units called *Integration Profiles* that highlight their capacity to address specific clinical needs.

Users familiar with these concepts can use integration statements as an aid to determine what level of integration this product supports with other systems, and what clinical and operational benefits such integration might provide.

IHE integration statement

PACSGEAR – MediaWriter IHE Integration Statement		Date
		15 May 2018
Vendor	Product Name	Version
PACSGear	Media Writer	5.1.9
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Portable Data for imaging	Portable Media Creator	
	Portable Media Importer	
Import Reconciliation Workflow	Importer	Unscheduled Import - No MWL Scheduled Import - Uses MWL
<u>Internet address for vendor's IHE information</u> : https://docs.hyland.com/		
Links to Standards Conformance Statements for the Implementation		
HL7		
DICOM	www.docs.hyland.com/PACSGEAR/MW/en_US/5.1/Print/MW_DCS.pdf	
Links to general information on IHE		
In North America: www.rsna.org/ihe	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j