



Media Writer

Technical Specifications

LX-DOC-MW5.1.9-TS-EN-REVC

Regulations and Compliance

Australian Sponsor: Emergo Australia Level 20 Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia	Brazil Sponsor: Media Writer Software Registro Anvisa nº: 80117589010 Responsável Técnico: Luiz Levy Cruz Martins CRF- 42415 IMPORTADOR: Emergo Brazil Import Importação e Distribuição de Produtos Médicos Hospitalares Ltda. Avenida Francisco Matarazzo, 1.752, Salas 502/503, Água Branca, São Paulo-SP, CEP - 05001-200 - Cidade: Brasília / UF: DF CNPJ: 04.967.408/0001-98 e-mail: brazilvigilance@emergogroup.com		
EU Sponsor: <table border="1"><tr><td>EC</td><td>REP</td></tr></table> EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands	EC	REP	
EC	REP		



Hyland LLC

4309 Hacienda Drive, Suite 500
Pleasanton, CA 94588 USA

Tel: 1-844-535-1404

Email: TS_PACSGEAR@hyland.com

© 2018 Hyland. Hyland and the Hyland logo are trademarks of Hyland LLC, registered in the United States and/or other countries. All other trademarks are the property of their respective owners.

Table of Contents

Regulations and Compliance	2
What is Media Writer?	4
Product requirements	4
Antivirus exclusion recommendation	4
Windows patches and updates	4

What is Media Writer?

Media Writer is a program that physicians can use to burn studies onto CDs/DVDs or USB flash drives, and upload studies to image exchange services. This makes studies available to other physicians through Open Image Exchange, and to patients who can view their studies on a PC. With the Media Writer System CD/DVD burner, physicians can also design and print labels on the CDs and DVDs.

Product requirements

Specification	Description
Hardware	Epson DiscProducer, DiscPublisher
Operating Systems	Windows 7 x32/x64 SP1, Windows 10
Memory	4 GB (8 GB for the BTO object compression feature)
Disk Space	120 GB

Antivirus exclusion recommendation

This application supports passive antivirus scans. You must ensure that passive scanning of the server conforms to the security policy at the installation site. Real-time scanning is allowed, but you must configure the antivirus program to exclude the following folders.

- GET UPDATED FOLDER PATHS
- GET UPDATED FOLDER PATHS
- GET UPDATED FOLDER PATHS
- GET UPDATED FOLDER PATHS
- GET UPDATED FOLDER PATHS

Windows patches and updates

You can install all standard Microsoft security patches and updates on the system unless Lexmark notifies you not to do so.