



# Media Writer

## Technical Specifications

DOC-MW5.2-TS-EN-REVA

## Regulations and Compliance

<b>Australian Sponsor:</b> Emergo Australia Level 20 Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia	<b>Brazil Sponsor:</b> Media Writer Software Registro Anvisa nº: 80117580229 Responsável Técnico: Luiz Levy Cruz Martins CRF- 42415 IMPORTADOR: Emergo Brazil Import Importação e Distribuição de Produtos Médicos Hospitalares Ltda. Avenida Francisco Matarazzo, 1.752, Salas 502/503, Água Branca, São Paulo-SP, CEP - 05001-200 CNPJ: 04.967.408/0001-98 e-mail: brazilvigilance@ul.com		
<b>EU Sponsor:</b> <table border="1"><tr><td>EC</td><td>REP</td></tr></table> <b>EMERGO EUROPE</b> Prinsessegracht 20 2514 AP The Hague The Netherlands	EC	REP	
EC	REP		



**Hyland LLC**

4309 Hacienda Drive, Suite 500  
Pleasanton, CA 94588 USA

**For support, contact:**

**Tel: 1-844-535-1404**

**Email: [Key-pacsgear@hyland.com](mailto:Key-pacsgear@hyland.com)**

**Web: <https://community.hyland.com>**

© 2019 Hyland. Hyland and the Hyland logo are trademarks of Hyland LLC, registered in the United States and/or other countries. All other trademarks are the property of their respective owners.

## Symbols used in this document



European Economic Area certification mark



ISO 7000:3500 - Electronic Instructions For Use

To indicate on product or product packaging that relevant information for use of the product is available in electronic form rather than, or in addition to, printed paper form.



ISO 7000:3082 - Manufacturer

To identify the manufacturer of a product.



ISO 7000: 0434B – Caution

To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.

# Table of Contents

Regulations and Compliance ..... 2

*Symbols used in this document* ..... 3

**What is Media Writer?** ..... 5

**Product requirements** ..... 5

    Antivirus exclusion recommendation ..... 5

    Windows patches and updates ..... 5

## What is Media Writer?

Media Writer is a program that physicians can use to burn studies onto CDs/DVDs. With the Media Writer System CD/DVD burner, physicians can also design and print labels on the CDs and DVDs.

## Product requirements

Specification	Description
Hardware	Epson DiscProducer, DiscPublisher
Operating Systems	Windows 7 x32/x64 SP1, Windows 10
Memory	4 GB (8 GB for the BTO object compression feature)
Disk Space	120 GB

### Antivirus exclusion recommendation

This application supports passive antivirus scans. You must ensure that passive scanning of the server conforms to the security policy at the installation site. Real-time scanning is allowed, but you must configure the antivirus program to exclude the following folders.

- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Buffer
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Config
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Images
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\History
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Jobs
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Logs
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Reports

### Windows patches and updates

You can install all standard Microsoft security patches and updates on the system unless Hyland notifies you not to do so.