



Media Writer  
Technical Specifications

DOC-MW6.1-TS-EN-REVB

# Documentation Notice

Information in this document is subject to change without notice. The software described in this document is furnished only under a separate license agreement and may only be used or copied according to the terms of such agreement. It is against the law to copy the software except as specifically allowed in the license agreement. This document or accompanying materials may contain certain information which is confidential information of Hyland Software, Inc. and its affiliates, and which may be subject to the confidentiality provisions agreed to by you.

Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright law, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the express written permission of Hyland Software, Inc. or one of its affiliates.

Hyland, HXP, OnBase, Alfresco, Nuxeo, and product names are registered and/or unregistered trademarks of Hyland Software, Inc. and its affiliates in the United States and other countries. All other trademarks, service marks, trade names and products of other companies are the property of their respective owners.

© 2023 Hyland Software, Inc. and its affiliates.

The information in this document may contain technology as defined by the Export Administration Regulations (EAR) and could be subject to the Export Control Laws of the U.S. Government including for the EAR and trade and economic sanctions maintained by the Office of Foreign Assets Control as well as the export controls laws of your entity's local jurisdiction. Transfer of such technology by any means to a foreign person, whether in the United States or abroad, could require export licensing or other approval from the U.S. Government and the export authority of your entity's jurisdiction. You are responsible for ensuring that you have any required approvals prior to export.

## Regulatory and compliance



**Manufacturer**  
**Hyland Software, Inc.**  
12919 Earhart Ave  
Auburn, CA 95602

### **Support**

Telephone: 1-844-535-1404

Web: <https://community.hyland.com>

## General warning and precautions



To ensure patient safety and adhere to HIPAA Privacy Rule involving patient health information, the Healthcare provider shall assign an unique Accession Number and Medical Record Number (MRN) to identify individual patient's health record affiliated with the patient to prevent erroneous medical data and duplicate identification with other patients.

## Symbols glossary



3082 - Manufacturer: Indicates the manufacturer.



2497 - Date of Manufacture: Indicates the date when the medical device was manufactured.



2493 - Catalog number: Indicates the manufacturer's catalog number so the medical device can be identified.



2498 - Serial number: Indicates the manufacturer's serial number so that a specific medical device can be identified.



1641 - Consult instructions for use: Indicates the need for the user to consult the instructions for use.



0434A - Caution: Indicates the need for the user to consult the instructions for use for important information such as warnings and cautions.



2610 - Patient number: Indicates a unique number associated with an individual patient.



Quantity: Indicates the number of units.



Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

# Table of Contents

**Documentation Notice** ..... **2**

**Regulatory and compliance** ..... **3**

**General warning and precautions** ..... **3**

**Symbols glossary** ..... **4**

**What is Media Writer?** ..... **6**

**Product requirements** ..... **6**

    Antivirus exclusion recommendation ..... 6

    Windows patches and updates ..... 6

## What is Media Writer?

Media Writer is a program intended to enable physicians and clinicians to burn studies onto CDs, DVDs, or USB flash drives. This makes studies available to other physicians and patients, who can view their studies on a PC. With the Media Writer System CD/DVD burner, physicians can also design and print labels on the CDs and DVDs. Media Writer is not for diagnostic use.

## Product requirements

Specification	Description
Hardware	Epson DiscProducer, DiscPublisher
Operating Systems	Windows 10
Memory	4 GB (8 GB for the BTO object compression feature)
Disk Space	120 GB

## Antivirus exclusion recommendation

This application supports passive antivirus scans. You must ensure that passive scanning of the server conforms to the security policy at the installation site. Real-time scanning is allowed, but you must configure the antivirus program to exclude the following folders.

- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Buffer
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Config
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Images
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\History
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Jobs
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Logs
- C:\Program Files (x86)\PACSGEAR\MediaWriter\Service\Reports

## Windows patches and updates

You can install all standard Microsoft security patches and updates on the system unless Hyland notifies you not to do so.