



ModLink Web Forms

Release Notes

LX-DOC-MLF2.0.0-RN-EN-REVB

Regulations and Compliance

<p>Australian Sponsor: Emergo Australia Level 20 Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia</p>	<p>Brazil Sponsor: ModLink Web Forms Software Registro Anvisa nº: 80117589010 Responsável Técnico: Júlia Zema Parente Pinto / CRF – DF 2508 IMPORTADOR: Emergo Brazil Import Importação e Distribuição de Produtos Médicos Hospitalares Ltda. Avenida Francisco Matarazzo, 1.752, Salas 502/503, Água Branca, São Paulo-SP, CEP - 05001-200 - Cidade: Brasília / UF: DF CNPJ: 04.967.408/0001-98 e-mail: brazilvigilance@emergogroup.com</p>		
<p>EU Sponsor:</p> <table border="1"><tr><td>EC</td><td>REP</td></tr></table> <p>EMERGO EUROPE Prinsessegracht 20 2514 AP The Hague The Netherlands</p>	EC	REP	
EC	REP		



Hyland LLC

4309 Hacienda Drive, Suite 500
Pleasanton, CA 94588 USA

Tel: 1-844-535-1404

Email: TS_PACSGEAR@hyland.com

© 2018 Hyland. Hyland and the Hyland logo are trademarks of Hyland LLC, registered in the United States and/or other countries. All other trademarks are the property of their respective owners.

Table of Contents

Regulations and Compliance	2
Overview	4
ModLink Web Forms 2.0.4 (March 2018)	4
Updates	4
ModLink Web Forms 2.0.0 (November 2016)	4
New features	4
Updates	4
ModLink Web Forms 1.0.0 (December 2015)	4
New features	4

Overview

This document describes the features of ModLink Web Forms. Medical staff can use the application to open customized web forms, and then gather information from dictation systems and DICOM Worklist, add to or change the information, and then send the changes back to the dictation system. They can also send an image of the form to a DICOM source.

ModLink Web Forms 2.0.4 (March 2018)

Updates

This release provides the following update.

- Resolved database errors by increasing the JSON string limit.

ModLink Web Forms 2.0.0 (November 2016)

New features

This release provides the following new features.

- You can now share forms between servers.
- You can now open and automatically fill in a form directly from patient demographics screens in PACS/EHR.
- Added advanced forms with new design options including background images, free sizing and placement of fields, and field-independent labels.

Updates

This release provides the following update.

- Resolved an issue regarding service memory usage.

ModLink Web Forms 1.0.0 (December 2015)

This release includes the following version.

- 1.0.0

New features

This release provides the following new features.

- Create and manage forms for viewing and modifying patient data on dictation systems
- Manage images for use in forms
- Populate forms with DICOM Worklist lookup results
- Manage users
- Manage destinations, query sources, devices, and other DICOM settings
- Customize configurations by department
- Administrative tools including licensing, and audit, server, and job logging