



ModLink

Technical Specifications

DOC-ML25-TS-EN-REVA

Documentation Notice

Information in this document is subject to change without notice. The software described in this document is furnished only under a separate license agreement and may only be used or copied according to the terms of such agreement. It is against the law to copy the software except as specifically allowed in the license agreement. This document or accompanying materials may contain certain information which is confidential information of Hyland Software, Inc. and its affiliates, and which may be subject to the confidentiality provisions agreed to by you.

Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright law, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the express written permission of Hyland Software, Inc. or one of its affiliates.

Hyland, HXP, OnBase, Alfresco, Nuxeo, and product names are registered and/or unregistered trademarks of Hyland Software, Inc. and its affiliates in the United States and other countries. All other trademarks, service marks, trade names and products of other companies are the property of their respective owners.

© 2025 Hyland Software, Inc. and its affiliates.

The information in this document may contain technology as defined by the Export Administration Regulations (EAR) and could be subject to the Export Control Laws of the U.S. Government including for the EAR and trade and economic sanctions maintained by the Office of Foreign Assets Control as well as the export controls laws of your entity's local jurisdiction. Transfer of such technology by any means to a foreign person, whether in the United States or abroad, could require export licensing or other approval from the U.S. Government and the export authority of your entity's jurisdiction. You are responsible for ensuring that you have any required approvals prior to export.

.

Regulatory and compliance



Manufacturer
HylandSoftware, Inc.
12919 Earhart Ave
Auburn, CA 95602 USA

Support

Telephone: 1-844-535-1404

Web: <https://community.hyland.com>

General warnings and precautions



To ensure patient safety and adhere to HIPAA Privacy Rule involving patient health information, the Healthcare provider shall assign an unique Accession Number and Medical Record Number (MRN) to identify individual patient's health record affiliated with the patient to prevent erroneous medical data and duplicate identification with other patients.

Symbols glossary



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

3082 - Manufacturer: Indicates the medical device manufacturer.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2497 - Date of Manufacture: Indicates the date when the medical device was manufactured.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2493 - Catalog number: Indicates the manufacturer's catalog number so the medical device can be identified.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2498 - Serial number: Indicates the manufacturer's serial number so that a specific medical device can be identified.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

1641 - Consult instructions for use: Indicates the need for the user to consult the instructions for use.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

0434A - Caution: Indicates the need for the user to consult the instructions for use for important information such as warnings and cautions.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2610 - Patient number: Indicates a unique number associated with an individual patient.



Quantity: Indicates the number of units.



Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

Table of Contents

| | |
|---|----------|
| Documentation Notice | 2 |
| Regulatory and compliance | 3 |
| General warnings and precautions | 3 |
| Symbols glossary | 4 |
| ModLink Technical Specifications | 6 |
| Setup program requirements | 6 |
| Operating systems | 6 |
| Hardware | 6 |
| Memory | 6 |
| Hard disk | 6 |
| Network interface | 6 |
| Virtual host requirements | 6 |
| Display | 7 |
| SQL | 7 |
| ModLink Web Forms | 7 |
| Antivirus exclusion recommendation | 7 |
| Windows patches and updates | 7 |

ModLink Technical Specifications

ModLink is a program that automatically populates clinical worksheets and forms for dictation systems with values from modalities, saving manual dictation time.

Setup program requirements

Setup program requirements by operating system are provided in the table below. The setup program prompts you and then installs certain prerequisites if not detected on the destination PC.

| OS of Host Computer | Requirement |
|---------------------|--|
| Windows Server 2019 | <ul style="list-style-type: none"> • IIS • .NET Extensibility 4.7 • IIS Hostable Web Core • Cumulative Update for .NET Framework 4.7.2 and 4.8 |
| Windows Server 2022 | <ul style="list-style-type: none"> • IIS • .NET Extensibility 4.7 • IIS Hostable Web Core • Cumulative Update for .NET Framework 4.7.2 and 4.8 |

Operating systems

- Windows Server 2019
- Windows Server 2022

Hardware

- Xeon, 4 cores, 2 GHz or later

Memory

- 8 GB or more

Hard disk

- 150 GB of free space

Network interface

- 100 MB/1 GB
- The server requires a static MAC address.

Virtual host requirements

The minimum virtual host requirements are provided in the following table.

| Specification | Description |
|--------------------|---------------------------------|
| Virtual processors | 8 vCPU, 1.8 GHz |
| Memory | 10 GB |
| Hard disk | 1 GB free space |
| Network interface | 100 MB/1 GB, static MAC address |

Display

- 1280 x 1024 32-bit color

SQL

- Microsoft SQL Server 2019 or later

Note:

The ModLink SQL Database needs approximately 20mb per 1000 studies sent through it.

ModLink Web Forms

If using ModLink WebForms, you must upgrade your PACSgear Core Server installation to version 5.3 or later.

Antivirus exclusion recommendation

This application supports passive antivirus scans. You must ensure that passive scanning of the server conforms to the security policy at the installation site. Real-time scanning is allowed, but you must configure the antivirus program to exclude the following folders. These folders represent the default install location. If you chose to install to a different location please adjust the exclusions to match your installation.

- C:\Program Files (x86)\Pacsgear\ModLink\DATA
- C:\Program Files (x86)\Pacsgear\ModLink\Jobs
- C:\Program Files (x86)\Pacsgear\ModLink\logs
- C:\Program Files (x86)\Pacsgear\ModLink\RawCap
- C:\Program Files (x86)\Pacsgear\ModLinkWeb\ModlinkWebService\MODLINKFILES
- C:\Program Files (x86)\Pacsgear\ModLinkWeb\ModlinkWebService\MODLINKFILES_DATA
- C:\Program Files (x86)\Pacsgear\ModLinkWeb\ModlinkWebService\MODLINK_DATA
- C:\Program Files (x86)\Pacsgear\ModLinkWeb\ModlinkWebService\JOBS
- C:\Program Files (x86)\Pacsgear\ModLinkWeb\ModlinkWebService\RawCap

Windows patches and updates

You can install all standard Microsoft security patches and updates on the system unless your representative notifies you not to do so.