



PACS Scan

PACS Scan, PACS Scan XDS PACS Scan Film

Technical Specifications

DOC-PS6.1-TS-EN-REVA

Documentation Notice

Information in this document is subject to change without notice. The software described in this document is furnished only under a separate license agreement and may only be used or copied according to the terms of such agreement. It is against the law to copy the software except as specifically allowed in the license agreement. This document or accompanying materials may contain certain information which is confidential information of Hyland Software, Inc. and its affiliates, and which may be subject to the confidentiality provisions agreed to by you.

Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright law, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the express written permission of Hyland Software, Inc. or one of its affiliates.

Hyland, HXP, OnBase, Alfresco, Nuxeo, and product names are registered and/or unregistered trademarks of Hyland Software, Inc. and its affiliates in the United States and other countries. All other trademarks, service marks, trade names and products of other companies are the property of their respective owners.

© 2024 Hyland Software, Inc. and its affiliates.

The information in this document may contain technology as defined by the Export Administration Regulations (EAR) and could be subject to the Export Control Laws of the U.S. Government including for the EAR and trade and economic sanctions maintained by the Office of Foreign Assets Control as well as the export controls laws of your entity's local jurisdiction. Transfer of such technology by any means to a foreign person, whether in the United States or abroad, could require export licensing or other approval from the U.S. Government and the export authority of your entity's jurisdiction. You are responsible for ensuring that you have any required approvals prior to export.

Regulatory and compliance



Manufacturer
Hyland Software, Inc.
12919 Earhart Ave.
Auburn, CA 95602 USA

Support
Telephone: 1-844-535-1404
Web: <https://community.hyland.com>

General warnings and precautions



To ensure patient safety and adhere to HIPAA Privacy Rule involving patient health information, the Healthcare provider shall assign an unique Accession Number and Medical Record Number (MRN) to identify individual patient's health record affiliated with the patient to prevent erroneous medical data and duplicate identification with other patients.

Symbols glossary



3082 - Manufacturer



2497 - Date of Manufacture



2493 - Catalog number



2498 - Serial number



1641 - Consult instructions for use: Indicates the need for the user to consult the instructions for use.



0434A - Caution: Indicates the need for the user to consult the instructions for use for important information such as warnings and cautions.



2610 - Patient number: Indicates a unique number associated with an individual patient.



Quantity: Indicates the number of units.



Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

Table of Contents

Documentation Notice	2
Regulatory and compliance	3
General warnings and precautions	3
Symbols glossary	4
About PACS Scan	6
Technical Specifications	6
PACS Scan requirements	6
Anti virus exclusion recommendation	6
Windows patches and updates	6

About PACS Scan

PACS Scan is a software program that captures patient images, video, and demographics from DICOM and non-DICOM sources (such as paper, film, CDs, and other electronic media), associates the information with existing patient records, and sends it to PACS/VNA.

Technical Specifications

PACS Scan requirements

The following list provides the minimum hardware and operating system requirements.

Specification	Requirement
Processor	Intel Core 2 Duo
Memory	4 GB 8 GB or more is recommended when importing images and studies 4 GB or larger
Hard disk	PACS Scan: 75 MB
Operating system	Windows 10 and 11
Display	1024 x 768, 32-bit color
Network interface	10/100/1000 Mbps
Network protocol	Standard network connection to DICOM devices and RIS system (if available).
Scanner	Standard TWAIN-compliant scanning device

Anti virus exclusion recommendation

This application supports passive anti virus scans. You must ensure that passive scanning of the server conforms to the security policy at the installation site. Real-time scanning is allowed, but you must configure the anti virus program to exclude the following folders.

- C:\Program Files\PACSGEAR\PacsSCAN
- C:\Users\Public\Documents\zvprt50\PacsSCAN

For proper functionality, open port 8182 in your anti-virus protection.

Windows patches and updates

You can install all standard Microsoft security patches and updates on the system unless your support representative notifies you not to do so.