

# PACS Scan Mobile

## Release Notes

Version: 1.1.x

Written by: Product Knowledge, R&D

Date: October 2015

LX-DOC-PSM-1.1.0-RN-EN-REVA

© 2015 Lexmark International Technology, S.A. All rights reserved.

Lexmark is a trademark of Lexmark International Technology, S.A., or its subsidiaries, registered in the U.S. and/or other countries. All other trademarks are the property of their respective owners. No part of this publication may be reproduced, stored, or transmitted in any form without the prior written permission of Lexmark.

# Table of Contents

- PACS Scan Mobile 1.1.x (October 2015)..... 4**
  - Enhancements..... 4
- PACS Scan Mobile 1.0.x (April 2015) ..... 4**
  - About PACS Scan Server..... 4
  - About PACS Scan Mobile Administration Console ..... 4
  - Features..... 4

## PACS Scan Mobile 1.1.x (October 2015)

The following enhancements and updates are included in this release.

### Enhancements

- PACS Scan Mobile is now available on iTunes.
- Administrators can now provision devices for users by setting the IP address and SSL options.
- You can now include a backslash in an MRN for patient entry and lookup.

## PACS Scan Mobile 1.0.x (April 2015)

PACS Scan Mobile is an app for mobile devices that physicians and authorized staff can use to capture images and videos easily. You can use the built-in camera of the mobile device to capture images or videos of your patient and then attach them to patient studies. When you are finished capturing images and videos, you must send them to a selected query source for storage, but you have the option of editing the images first.

### About PACS Scan Server

PACS Scan Server is a desktop program that controls the flow of data between DICOM or related systems such as PACS, archives, Worklist servers, modalities, EHR/EMR, and RIS, and any clients that connect to it such as the PACS Scan Mobile app.

### About PACS Scan Mobile Administration Console

PACS Scan Mobile Administration Console (PSM Console) is a browser-based interface to PACS Scan Server that administrators can use to configure server settings, and settings for all of the client mobile devices that connect to the server. These include settings for licensing, DICOM querying and sending, users, and logging.

### Features

PACS Scan Mobile offers the following features.

- Query and display patient demographics from multiple sources, or create new patient data.
- Capture clinical images and videos, attach to patient studies, and send to DICOM archives.
- Edit demographics and media.
- Manage mobile device configuration and user access from a web console.
- Create department-specific configurations.
- Make both DICOM and Worklist queries, and send HL7® messages.
- View DICOM Trace data, server logs, and audit logs.