



# PACS Scan Mobile

## Technical Specifications


LX-DOC-PSM3.0-TS-EN-REVA

## Regulations and Compliance

**Australian Sponsor:**

Emergo Australia  
Level 20 Tower II  
Darling Park 201 Sussex Street  
Sydney, NSW 2000 Australia

**EU Sponsor:**



**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



**Manufacturer**

**Hyland LLC**

4309 Hacienda Drive, Suite 500  
Pleasanton, CA 94588 USA

**For support, contact:**

**Tel: 1-844-535-1404**

**Web: <https://community.hyland.com>**

© 2019 Hyland. Hyland and the Hyland logo are trademarks of Hyland LLC, registered in the United States and/or other countries. All other trademarks are the property of their respective owners.

## Symbols used in this document



CE marking of conformity

The device is in conformity with the applicable requirements set out in this EU Medical Device Regulation and other applicable Union harmonization legislation.



ISO 7000:3500 - Electronic Instructions For Use

To indicate on product or product packaging that relevant information for use of the product is available in electronic form rather than, or in addition to, printed paper form.



ISO 7000:3082 - Manufacturer

To identify the manufacturer of a product.



ISO 7000: 0434B – Caution

To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.

# Table of Contents

<b>Regulations and Compliance</b> .....	<b>2</b>
Symbols used in this document .....	3
<b>What is PACS Scan Mobile?</b> .....	<b>5</b>
<b>Technical specifications</b> .....	<b>5</b>
PACS Scan Mobile app requirements .....	5
PACSGear Core Server requirements .....	5

## What is PACS Scan Mobile?

The PACS Scan Mobile client is an iOS application for mobile devices that allows you to easily capture images and videos and then send them to patient studies on DICOM archives. With proper authorization, you can query and retrieve a patient study from a worklist or create a new patient record. You can take photos and videos. You can manipulate photos, for example, crop or rotate, and also annotate them prior to sending to other locations. You cannot alter videos prior to sending. You can delete unwanted images and videos prior to sending.

## Technical specifications

### PACS Scan Mobile app requirements

The minimum hardware and operating system requirements are provided in the following table.

Specification	Description
Hardware	iPhone <ul style="list-style-type: none"> <li>• X</li> <li>• 8 Plus</li> <li>• 8</li> <li>• SE</li> <li>• 7 Plus</li> <li>• 7</li> <li>• 6s Plus</li> <li>• 6s</li> </ul> iPad <ul style="list-style-type: none"> <li>• Pro (10.5", 12.9" and older models of 9.7"</li> <li>• 2017</li> </ul>
Operating Systems	iOS 11 and later
Web Application Servers	IIS 7
Memory	64 MB
Internet Browser	Safari

### PACSGear Core Server requirements

PACS Scan Mobile is included with PACSGear Core Server. Your server must meet the requirements described in the PACSGear Core Server technical specifications.