



PACS Scan Mobile
Technical Specifications

DOC-PSM3.1-TS-EN-REVA

Copyright

Information in this document is subject to change without notice. The software described in this document is furnished only under a separate license agreement and may be used or copied only according to the terms of such agreement. It is against the law to copy the software except as specifically allowed in the license agreement. This document or accompanying materials contains certain information which is confidential information of Hyland Software, Inc. and its affiliates, and which is subject to the confidentiality provisions agreed to by you.

All data, names, and formats used in this document's examples are fictitious unless noted otherwise. Complying with all applicable copyright laws is the responsibility of the user. Without limiting the rights under copyright law, no part of this document may be reproduced, stored in or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the express written permission of Hyland Software, Inc. or one of its affiliates.

Hyland®, Hyland Software®, Hyland Healthcare, and Hyland product names, are registered and/or unregistered trademarks of Hyland Software, Inc. and its affiliates in the United States and other countries. All other trademarks, service marks, trade names and products of other companies are the property of their respective owners.

© 2020 Hyland Software, Inc. and its affiliates. All rights reserved.

Regulatory and compliance

Australian Sponsor:

Emergo Australia
Level 20 Tower II
Darling Park 201 Sussex Street
Sydney, NSW 2000 Australia

EU Representative:**EMERGO EUROPE**

Prinsessegracht 20
2514 AP The Hague
The Netherlands

**Manufacturer****Hyland LLC**

4309 Hacienda Drive, Suite 500
Pleasanton, CA 94588 USA

Support

Telephone: 1-844-535-1404

Web: <https://community.hyland.com>

General warning and precautions



To ensure patient safety and adhere to HIPAA Privacy Rule involving patient health information, the Healthcare provider shall assign a unique Accession Number and Medical Record Number (MRN) to identify individual patient's health record affiliated with the patient to prevent erroneous medical data and duplicate identification with other patients.

Symbols glossary



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

3082 - Manufacturer: Indicates the medical device manufacturer.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

Authorized representative in European Community: Indicates the authorized representative in the European union. Symbol is accompanied by the name and address of the authorized representative adjacent to the symbol.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2497 - Date of Manufacture: Indicates the date when the medical device was manufactured.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2493 - Catalog number: Indicates the manufacturer's catalog number so the medical device can be identified.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2498 - Serial number: Indicates the manufacturer's serial number so that a specific medical device can be identified.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

1641 - Consult instructions for use: Indicates the need for the user to consult the instructions for use.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

0434A - Caution: Indicates the need for the user to consult the instructions for use for important information such as warnings and cautions.



ISO - 15223-1 - Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied.

2610 - Patient number: Indicates a unique number associated with an individual patient.



Medical Devices Directive 93/42/EEC

CE Marking: Indicates that a product complies with applicable European Union regulations.



Quantity: Indicates the number of units.

Table of Contents

Copyright	2
Regulatory and compliance	3
General warning and precautions	3
Symbols glossary	4
What is PACS Scan Mobile?	6
Technical specifications	6
PACS Scan Mobile app requirements	6
PACSGear Enterprise Server requirements	6

What is PACS Scan Mobile?

The PACS Scan Mobile client is an iOS application for mobile devices that allows you to easily capture images and videos and then send them to patient studies on DICOM archives. With proper authorization, you can query and retrieve a patient study from a worklist or create a new patient record. You can take photos and videos. You can manipulate photos, for example, crop or rotate, and also annotate them prior to sending to other locations. You cannot alter videos prior to sending. You can delete unwanted images and videos prior to sending.

Technical specifications

PACS Scan Mobile app requirements

The minimum hardware and operating system requirements are provided in the following table.

Specification	Description
Hardware	iPhone <ul style="list-style-type: none"> • X • 8 Plus • 8 • SE • 7 Plus • 7 • 6s Plus • 6s iPad <ul style="list-style-type: none"> • Pro (10.5", 12.9" and older models of 9.7" • 2017
Operating Systems	iOS 11 and later
Web Application Servers	IIS 7
Memory	64 MB
Internet Browser	Safari

PACSGear Enterprise Server requirements

PACS Scan Mobile is included with PACSGear Enterprise Server. Your server must meet the requirements described in the PACSGear Enterprise Server technical specifications.