

Perceptive Audit Record Repository Manager

Release Notes

Version: 1.0.x

Written by: Product Knowledge, R&D

Date: Friday, November 04, 2016

©2016 Lexmark. All rights reserved.

Lexmark is a trademark of Lexmark International Inc., registered in the U.S. and/or other countries. All other trademarks are the property of their respective owners. No part of this publication may be reproduced, stored, or transmitted in any form without the prior written permission of Lexmark.

Table of Contents

Perceptive Audit Record Repository Manager Release (August, 2016) 4

Perceptive Audit Record Repository Manager Release (August, 2016)

This release includes the following version.

- Perceptive Audit Record Repository Manager 1.0.x

The Perceptive Audit Record Repository Manager, version 1.0.x, a central audit server repository, provides secure (TLS TCP) and non-secure (UDP and TCP) transmission and auditing options for different healthcare application systems.

Audit Record Repository (ARR) is a system unit that receives and collects audit records from multiple systems and helps to meet the requirements for patient information privacy as prescribed by the Health Insurance Portability and Accountability Act (HIPAA) regulation. Exchanges and access to this information that is outside the scope of regular treatment, planning, or healthcare operations should be registered in an audit trail that can be regularly inspected. Integrating the Healthcare Enterprise (IHE) specifies the exchange protocol between individual devices that collect the audit trail information using a central audit repository to replace multiple audit trails for each individual system.

The IHE initiative defines a set of profiles, actors, and transactions that create common scenarios for particular workflow processes. The Integration Profiles depict security as a fundamental requirement of the framework. Specifically, the Audit Trail and Node Authentication (ATNA) profile defines standards-based mechanisms for securely transmitting and storing audit records in a central repository. The data structure defined by the profile provides a number of record types that capture different audit events.

The Perceptive Audit Record Repository Manager, which conforms to the IHE standards, provides a mechanism for efficient searching, viewing, and deletion of ATNA messages. It can receive messages in both RFC-3881 and the new DICOM schemas.